

# SCHEIN ORTHOPÄDIE SERVICE



**WeXhelp**

INCREASING QUALITY OF LIFE  
AND WELLBEING OF PEOPLE!

# Index

## Children's shoes

<b>My generation</b>	Orthotic shoes	4–19
	Support shoes	20–25
	Multi-support shoes	26–27
	Anti-varus/Neutral shoes	28–33

## Comfort shoes

<b>Lucro</b> Special shoes for people with diabetes or rheumatism	Classic	34–37
	Kinetic	38–39
	Sportic	40–41
	Ergonic	42–43
<b>Actiflex</b>	The feel-good shoes	44–47
<b>Ortho lady &amp; gent</b>	Comfort shoes	48–51

## Post-operative care

<b>Post op shoes</b>	Shoes for post-surgery care	52–55
----------------------	-----------------------------	-------

## Insoles

<b>Novaped</b>	Performance	56–57
<b>Novaped</b>	Carbon	58–60
<b>Novaped</b>	Elements	61
<b>Novaped</b>	Heel plateau	62–67
<b>Novaped</b>	Insoles	68–69
<b>Novaped</b>	Moulding sets	70–71
<b>Novaped</b>	Sports	72–77
<b>Novaped</b>	Soft	78–85
<b>Novaped</b>	Soft–Ideal for diabetics	86–87
<b>Novaped</b>	Blanks–Ideal for diabetics	88–89

## Schein materials

<b>Schein materials</b>	Schein padding and bedding materials	90–99
<b>Schein materials</b>	Schein/ARU® combination plates	100–103
<b>Schein materials</b>	Schein thermoplastics and cork	104–105

## Space for therapy

Many people often find wearing a leg orthosis restricts their freedom of movement. Our orthotic shoes are therefore designed to ensure an optimal fit for the foot, orthosis and therapeutic shoe. The sturdy heel cap in combination with the raised and flexible shaft edge ensures a good fit and high stability.

But an ideal fit is not all that matters. In everyday life in particular, it is extremely important that our orthotic shoes are easy to put on and take off.



With that in mind, many of our models have a tab to simplify putting on and removing them. But that's not all: the extended tongue opening in many models makes it possible to open the entire shoe very wide.



In this Tutorials answer we all questions concerning our orthopaedic shoe range.

Enjoy the videos!

If further adjustments are required for the orthosis, they can be made individually by an orthopaedic technician. Specific orthopaedic measures and additions to ensure the feet are optimally cared for can be implemented on all orthotic shoes.



Available in the widths TN5, TN8, TN9 and TN12.

The widths and sizes of shoes can be combined in cases of one-sided treatment.

Diagram: symbolic representation of insole shape and last volume.

The models shown in this brochure are a choice of our collection.

In order to see Schein's complete shoe range please go to: <https://mygeneration-schein.de/en/>





Tuffy

Art. no.	Width	Size
371090	TN8	24-42



You can also design  
this model yourself in the  
configurator

my design

www.mydesign-schein.de



Domino

Art. no.	Width	Size
371050	TN8	27-42



Hero

Art. no.	Width	Size
361080	TN5	18-45
371080	TN8	18-45



Chen

Art. no.	Width	Size
367000	TN5	20-45
377000	TN8	20-45
377005	TN12	20-45



Samara

Art. no.	Width	Size
367060	TN5	21-32
377060	TN8	21-32



Paula

Art. no.	Width	Size
367030	TN5	20-38
377030	TN8	20-38



Nell

Art. no.	Width	Size
367070	TN5	24-35
377070	TN8	24-35

Our orthopaedic shoes are registered at NAV Norway for refund.  
For verification please go to: <https://lovdata.no/nav/rundskriv/v2-10-07i>

## Orthotic shoes



### Crossi

Art. no.	Width	Size
361950	TN5	21-32
371950	TN8	21-32



### Kletty

Art. no.	Width	Size
363950	TN5	21-35
373950	TN8	21-35



### Mareike

Art. no.	Width	Size
363960	TN5	21-32
373960	TN8	21-32



### Carbon

Art. no.	Width	Size
361900	TN5	21-45
371900	TN8	21-45



### Ocean

Art. no.	Width	Size
362550	TN5	18-35
372550	TN8	18-35



### Cora

Art. no.	Width	Size
364630	TN5	21-35
374630	TN8	21-35



### Pecore

Art. no.	Width	Size
364680	TN5	24-39
374680	TN8	24-39

Our orthopaedic shoes are registered at NAV Norway for refund.  
For verification please go to: <https://lovdata.no/nav/rundskriv/v2-10-07i>



## Navaro

Art. no.	Width	Size
364200	TN5	30–47
374200	TN8	30–47



## Rio

Art. no.	Width	Size
364210	TN5	24–39
374210	TN8	24–39



## Bailey

Art. no.	Width	Size
364260	TN5	30–42
374260	TN8	30–42



## Fynn

Art. no.	Width	Size
362810	TN5	30–47
370810	TN8	30–47



## Paco

Art. no.	Width	Size
366720	TN5	20–47
376720	TN8	20–47
376725	TN9	20–47
378725	TN12	20–47




## Fox

Art. no.	Width	Size
366780	TN5	21–47
376780	TN8	21–47



## Van

Art. no.	Width	Size
368700	TN5	36–50
378700	TN8	36–50



Our orthopaedic shoes are registered at NAV Norway for refund.  
For verification please go to: <https://lovdata.no/nav/rundskriv/v2-10-07i>

## Orthotic shoes



### Blue

Art. no.	Width	Size
363320	TN5	21-42
373320	TN8	21-42



### Buddy

Art. no.	Width	Size
363750	TN5	26-42
373750	TN8	26-42



### Maya

Art. no.	Width	Size
363310	TN5	21-38
373310	TN8	21-38



### Blacky

Art. no.	Width	Size
363030	TN5	30-47
373030	TN8	30-47



### Gym

Art. no.	Width	Size
363010	TN5	26-47
373010	TN8	26-47



### Theo

Art. no.	Width	Size
369600	TN5	30-45
379600	TN8	30-45

Our orthopaedic shoes are registered at NAV Norway for refund.  
For verification please go to: <https://lovdata.no/nav/rundskriv/v2-10-07i>



## Alaska

Art. no.	Width	Size
369420	TN5	27-39
379420	TN8	27-39



## Conny

Art. no.	Width	Size
369430	TN5	27-35
379430	TN8	27-35



## Orthotic shoes



### Hoover

Art. no.	Width	Size
369580	TN5	24-42
379580	TN8	24-42



### Rover

Art. no.	Width	Size
369500	TN5	27-39
379500	TN8	27-39



### Noel

Art. no.	Width	Size
369780	TN5	24-39
379780	TN8	24-39



### Adrian

Art. no.	Width	Size
369700	TN5	24-39
379700	TN8	24-39



### Nanaimo

Art. no.	Width	Size
368220	TN5	24-42
378220	TN8	24-42



### Alina

Art. no.	Width	Size
368280	TN5	24-39
378280	TN8	24-39



### Trivor

Art. no.	Width	Size
368200	TN5	36-45
378200	TN8	36-45

Our orthopaedic shoes are registered at NAV Norway for refund.  
For verification please go to: <https://lovdata.no/nav/rundskriv/v2-10-07i>



## Dori

Art. no.	Width	Size
364700	TN5	27-42
374700	TN8	27-42



## Pongo

Art. no.	Width	Size
364720	TN5	21-35
374720	TN8	21-35




## Lucky

Art. no.	Width	Size
364780	TN5	21-35
374780	TN8	21-35



## Penny

Art. no.	Width	Size
364730	TN5	21-32
374730	TN8	21-32



## Orthotic shoes



### Viva

Art. no.	Width	Size
364370	TN5	33–45
374370	TN8	33–45



### Rose

Art. no.	Width	Size
364380	TN5	21–35
374380	TN8	21–35



### Leon

Art. no.	Width	Size
364310	TN5	27–47
374310	TN8	27–47



### Dom

Art. no.	Width	Size
367180	TN5	27–39
377180	TN8	27–39



### Perdita

Art. no.	Width	Size
367130	TN5	27–35
377130	TN8	27–35



### Tuvato

Art. no.	Width	Size
367100	TN5	30–42
377100	TN8	30–42

Our orthopaedic shoes are registered at NAV Norway for refund.  
For verification please go to: <https://lovdata.no/nav/rundskriv/v2-10-07i>





## Milow

Art. no.	Width	Size
369200	TN5	21-47
379200	TN8	21-47



## Jewel

Art. no.	Width	Size
369280	TN5	24-39
379280	TN8	24-39



You can also design  
this model yourself in the  
configurator

**my design**

[www.mydesign-schein.de](http://www.mydesign-schein.de)



## Orthotic shoes



### Missy

Art. no.	Width	Size
369830	TN5	30–39
379830	TN8	30–39



### Hammy

Art. no.	Width	Size
369800	TN5	30–47
379800	TN8	30–47



### Silver

Art. no.	Width	Size
369880	TN5	30–42
379880	TN8	30–42



### Curly

Art. no.	Width	Size
364110	TN5	20–42
374110	TN8	20–42



### Jacko

Art. no.	Width	Size
364180	TN5	23–38
374180	TN8	23–38



You can also design  
this model yourself in the  
configurator

my design

[www.mydesign-schein.de](http://www.mydesign-schein.de)



### Maris

Art. no.	Width	Size
364140	TN5	23–45
374140	TN8	23–45



### Julian

Art. no.	Width	Size
364170	TN5	20–42
374170	TN8	20–42
374175	TN12	20–42



## Ole

Art. no.	Width	Size
364920	TN5	21–32
374920	TN8	21–32



## River

Art. no.	Width	Size
364980	TN5	30–47
374980	TN8	30–47



## Dani

Art. no.	Width	Size
363460	TN5	24–42
373460	TN8	24–42



## Vegas

Art. no.	Width	Size
363420	TN5	24–45
373420	TN8	24–45



## Phoenix

Art. no.	Width	Size
363480	TN5	18–50
373480	TN8	18–50
373485	TN9	18–50
378485	TN12	18–50



## Stella

Art. no.	Width	Size
363450	TN5	21–38
373450	TN8	21–38
373455	TN9	21–38
378455	TN12	21–38

## Orthotic shoes



Micki

Art. no.	Width	Size
363530	TN5	36-47



Robin

Art. no.	Width	Size
363510	TN5	36-47
373510	TN8	36-47



Matterhorn

Art. no.	Width	Size
363620	TN5	30-47





### Nevis

Art. no.	Width	Size
362600	TN5	20-47
372600	TN8	20-47



Valid for all models  
on this page



### Berni

Art. no.	Width	Size
372510	TN8	26-38



### Felix

Art. no.	Width	Size
362500	TN5	26-42
372500	TN8	26-42

## The perfect twist!

Originally developed for sport, this quick-lacing system for orthotics wearers with limited motor skills makes putting on and removing shoes considerably easier. The loosening and tightening is performed with the turning wheel. When the wheel is turned, all the laces are uniformly adjusted in a manner adapted to the patient's foot and set with accuracy down to the millimetre. This fastening ensures a stable hold and the highest level of comfort when wearing the shoe.

The sturdy design of the atop fastening system supports the philosophy of our orthotic shoes:

- Very good hold
- Fine adjustment
- One-handed operation
- Low weight
- Available as a front or rear fastening



## Secure hold

We offer two versions of the atop fastening system:

- With turning wheel on the rear cap
- With turning wheel on the middle of the shoe

Both lacing systems provide excellent hold by uniformly tightening the entire surface.

Compared to conventional laces, which can absorb water or dirt, the low-friction plastic covers wrapped around the fine stainless-steel wires maintain their featherlight weight.

Even under heavy strain, the coil and wire form a long-lasting combination.

Once affixed, the lacing stays in place—the orthotics or the foot are held securely.



### To loosen:

simply turn the wheel to the left and the laces are loosened.



### To tighten:

turn the wheel to the right and the laces are fixed tightly.

The models shown in this brochure are a choice of our collection. In order to see Schein's complete shoe range please follow: <https://mygeneration-schein.de/en/>



### Minerva

Art. no.	Width	Size
362940	TN5	24–32
372940	TN8	24–32



### Juno

Art. no.	Width	Size
366250	TN5	24–39
376250	TN8	24–39



### Apollo

Art. no.	Width	Size
360040	TN5	24–47
370040	TN8	24–47



### Äther

Art. no.	Width	Size
362970	TN5	24–47
372970	TN8	24–47

Our orthopaedic shoes are registered at NAV Norway for refund.  
For verification please go to: <https://lovdata.no/nav/rundskriv/v2-10-07i>



**Your trusty companion—**  
with flexible and individual width adjustment

Our aim is to form a strong team with our customers, doctors, physiotherapists, health insurance companies and medical services and to develop the best assistive technology in cooperation with our customers.

Our support shoes meet all necessary requirements for recognition by health insurance companies as therapeutic aids. The widths and sizes of shoes can easily be combined in cases of one-sided treatment.



Important criteria for the quality of shoes are:

- Soft materials on the shaft ends to avoid pressure points and ensure high levels of comfort when wearing
- High level of heel stabilisation with the foot free to roll at the same time
- Hook-and Loop fasteners and zippers to make the shoes easier to put on and remove



In this Tutorial answer we all questions concerning our orthopaedic shoe range.

**Enjoy the videos!**



Available in the sizes TN2, TN4 and TN6.

The widths and sizes of shoes can be combined in cases of one-sided treatment.

Diagram: symbolic representation of insole shape and last volume.

The models shown in this brochure are a choice of our collection.

In order to see Schein's complete shoe range please go to: <https://mygeneration-schein.de/en/>





Joshua

Art. no.	Width	Size
356820	TN4	27–35



Hero

Art. no.	Width	Size
351085	TN2	23–42
351080	TN4	23–42



Sunny

Art. no.	Width	Size
350450	TN4	23–35



Tuffy

Art. no.	Width	Size
351090	TN4	30–42



Christin

Art. no.	Width	Size
351730	TN2	23–35



Nicky

Art. no.	Width	Size
351700	TN2	23–35



Dori

Art. no.	Width	Size
354700	TN2	36–42



Our orthopaedic shoes are registered at NAV Norway for refund.  
For verification please go to: <https://lovdata.no/nav/rundskriv/v2-10-07i>

## Support shoes



Leon

Art. no.	Width	Size
354310	TN2	24-45



Piri

Art. no.	Width	Size
354360	TN2	30-42



Kaleb

Art. no.	Width	Size
350100	TN4	30-42



Jarle

Art. no.	Width	Size
359420	TN2	23-47

Our orthopaedic shoes are registered at NAV Norway for refund.  
For verification please go to: <https://lovdata.no/nav/rundskriv/v2-10-07i>



Dom

Art. no.	Width	Size
357180	TN2	33-39



Tuvato

Art. no.	Width	Size
357100	TN2	36-42



Silver

Art. no.	Width	Size
359880	TN2	36-42



## Support shoes



River

Art. no.	Width	Size
354980	TN2	30–47



my design

www.mydesign-schein.de



Unique

Art. no.	Width	Size
352680	TN2	27–42



Julian

Art. no.	Width	Size
354170	TN2	23–42



Stella

Art. no.	Width	Size
353450	TN2	27–39
353455	TN4	27–39



Phoenix

Art. no.	Width	Size
354480	TN2	23–45
353480	TN4	23–45
353485	TN6	23–42



Dani

Art. no.	Width	Size
353460	TN2	23–42
353465	TN4	23–42



Vegas

Art. no.	Width	Size
363420	TN5	24–45
373420	TN8	24–45

Our orthopaedic shoes are registered at NAV Norway for refund.  
For verification please go to: <https://lovdata.no/nav/rundskriv/v2-10-07i>



## Matterhorn

Art. no.	Width	Size
353620	TN2	30-47
353625	TN4	30-47



## Nevis

Art. no.	Width	Size
352600	TN2	23-47



**Higher shaft–**  
with a stable arthrodesis cap

These models can be further reinforced by two thermoplastically mouldable stabilising rods on the inner and outer side. The stiffening elements are available in two strengths: white = strong, skin-coloured = medium strength.

Four white rods are included.

This allows for even more lateral stability and rotational stability in the ankle joint area, without com-promising the rolling of the foot. The integrated rolling section on the ball of the foot facilitates step initiation.

As is the case with our other therapeutic shoes, individual shoe adjustments can also be made to these models. The applications for the multi-support shoe include, among other things, the follow-up treatment of conservatively or surgically treated capsular ligament lesions in the ankle, capsular ligament insufficiency or Achilles tendon injuries. It is also suitable for immobilisation of the tarsal and metatarsal areas, as well as follow-up treatment of fractures and arthrodesis.

Optional: Lateral rods provide additional stiffening and stability

High over-the-ankle cap to stabilise the ankle joint

Practical fastening and width adjustment with Hook-and-Loop fastener straps or laces

Shaft and strap padding



The models shown in this brochure are a choice of our collection. In order to see Schein's complete shoe range please go to: <https://mygeneration-schein.de/en/>



### Nevada

Art. no.	Width	Size
359000	ST3	36-47



### Montana

Art. no.	Width	Size
359010	ST3	36-47



### Oregon

Art. no.	Width	Size
359500	ST3	36-47





## Anti-varus shoes– for secure correction right from the start

At first, many children find it difficult to adjust to walking in anti-varus shoes. The stiffness of the shoe, which is a must from a therapeutic standpoint, is unaccustomed and makes walking more difficult because the foot is forced into the corrected position. However, experience has shown that children quickly become accustomed to the new feeling.

It is important to continuously monitor the treatment progress in consultation with your doctor. The earlier your child receives the right care, the easier it will be to rectify the problem. Where each foot requires different care, shoes of differing widths and sizes can be combined.


In cases of slight misalignment, therapy with anti-varus shoes can be supported by the following:

- Regular foot massage
- Suitable exercises
- Physiotherapy



In this Tutorial answer we all questions concerning our orthopaedic shoe range.

Enjoy the videos!



The widths and sizes of shoes can be combined in cases of one-sided treatment.

The models shown in this brochure are a choice of our collection.

In order to see Schein's complete shoe range please go to: <https://mygeneration-schein.de/en/>





Ben

Art. no.	Width	Size
041750	AV	18-42
041755	NS	18-42



Anja

Art. no.	Width	Size
041790	AV	18-28
041795	NS	18-28



Hero

Art. no.	Width	Size
041080	AV	18-42
041085	NS	18-42



Emilia

Art. no.	Width	Size
343040	AV	21-30
343045	NS	21-30



Lupi

Art. no.	Width	Size
041010	AV	18-42
041015	NS	18-42



Amy

Art. no.	Width	Size
041030	AV	18-32
041035	NS	18-32



Luik

Art. no.	Width	Size
341980	AV	21-35
341985	NS	21-35



Mareike

Art. no.	Width	Size
041620	AV	21-32
041625	NS	21-32

Our orthopaedic shoes are registered at NAV Norway for refund.  
For verification please go to: <https://lovdata.no/nav/rundskriv/v2-10-07i>

## Anti-varus/Neutral shoes



Pepper

Art. no.	Width	Size
347160	AV	21-29
347165	NS	21-29



Dom

Art. no.	Width	Size
347180	AV	21-29
347185	NS	21-29



Leon

Art. no.	Width	Size
344310	AV	21-42
344315	NS	21-42



Rose

Art. no.	Width	Size
344380	AV	18-29
344385	NS	18-29



Missy

Art. no.	Width	Size
349830	AV	21-29
349835	NS	21-29



Hammy

Art. no.	Width	Size
349800	AV	21-32
349805	NS	21-32

Our orthopaedic shoes are registered at NAV Norway for refund.  
For verification please go to: <https://lovdata.no/nav/rundskriv/v2-10-07i>



## Rolly

Art. no.	Width	Size
354710	AV	21-29
354715	NS	21-29

## Lucky

Art. no.	Width	Size
354780	AV	21-29
354785	NS	21-29



## Kerry

Art. no.	Width	Size
344890	AV	20-38
344895	NS	20-38



## Anti-varus/Neutral shoes



Ole

Art. no.	Width	Size
344920	AV	21–29
344925	NS	21–29



Luana

Art. no.	Width	Size
344910	AV	21–29
344915	NS	21–29



Julian

Art. no.	Width	Size
344170	AV	20–38
344175	NS	20–38



Momo

Art. no.	Width	Size
348120	AV	21–32
348125	NS	21–32



Vegas

Art. no.	Width	Size
343420	AV	24–39
343425	NS	24–39



Phoenix

Art. no.	Width	Size
042560	AV	18–42
042565	NS	18–42



Stella

Art. no.	Width	Size
343450	AV	18–32
343455	NS	18–32

Our orthopaedic shoes are registered at NAV Norway for refund.  
For verification please go to: <https://lovdata.no/nav/rundskriv/v2-10-07i>



### Nevis

Art. no.	Width	Size
342600	AV	20–42
342605	NS	20–42



### Cloud

Art. no.	Width	Size
042030	AV	18–42
042035	NS	18–42




### Marty

Art. no.	Width	Size
042050	AV	20–42
042055	NS	20–42



### Maxi

Art. no.	Width	Size
042040	AV	18–38
042045	NS	18–38



## The classic one

The LucRo classic can protect against recurring ulcers in diabetic foot syndrome with almost 80% certainty.

This was clinically proven in a 42-month study.

### Women



**Monika**

Art. no.	Width ex stock	Delivery time four weeks	Size
22423	S, M, W	Upper	36-43



**Monika**

Art. no.	Width ex stock	Delivery time four weeks	Size
22320	S, M, W	Upper	34-45



**Anita**

Art. no.	Width ex stock	Delivery time four weeks	Size
22409	M	S, W, Upper	36-43



**Marisa**

Art. no.	Width ex stock	Delivery time four weeks	Size
22065	M, W	S, +10, Upper	36-43



**Sylvia**

Art. no.	Width ex stock	Delivery time four weeks	Size
22040	S, M, W, +10, Upper		34-45



**Sylvia**

Art. no.	Width ex stock	Delivery time four weeks	Size
22445	M	S, W, +10, Upper	36-43

The models shown in this brochure are a choice of our collection.  
In order to see Schein's complete shoe range please go to: <https://www.lucro.de/en/>



Ilona

Art. no.	Width ex stock	Delivery time four weeks	Size
22362	M, W, +10	S, Upper	34-45



Ilona

Art. no.	Width ex stock	Delivery time four weeks	Size
22363	M	S, W, +10, Upper	36-43



Yvonne

Art. no.	Width ex stock	Delivery time four weeks	Size
22300	M, W	S, +10, Upper	36-43



Andrea

Art. no.	Width ex stock	Delivery time four weeks	Size
22190	M, W	S, +10, Upper	36-43



Doris

Art. no.	Width ex stock	Delivery time four weeks	Size
22150	M, W	S, +10, Upper	36-43



Hildegard

Art. no.	Width ex stock	Delivery time four weeks	Size
22394	S, M, W, +10	Upper	36-43

The models shown in this brochure are a choice of our collection.  
 In order to see Schein's complete shoe range please go to: <https://www.lucro.de/en/>



# classic

## Men



Eric

Art. no.	Width ex stock	Delivery time four weeks	Size
21276	S, M, W	Upper	40-47



Achim

Art. no.	Width ex stock	Delivery time four weeks	Size
21190	M, W, Upper	S	40-47



Peter

Art. no.	Width ex stock	Delivery time four weeks	Size
21200	M, W, Upper	S, +10	40-47



Oliver

Art. no.	Width ex stock	Delivery time four weeks	Size
21095	M, W, +10	S, Upper	40-47



Volker

Art. no.	Width ex stock	Delivery time four weeks	Size
21330	M, W, +10	S, Upper	37-50



Walter

Art. no.	Width ex stock	Delivery time four weeks	Size
21370	M, W, +10	S, Upper	37-50



Sven

Art. no.	Width ex stock	Delivery time four weeks	Size
21050	S, M, W, +10, Upper		37-50



Sven

Art. no.	Width ex stock	Delivery time four weeks	Size
21053	S, M, W, +10, Upper		37-50

The models shown in this brochure are a choice of our collection.  
In order to see Schein's complete shoe range please go to: <https://www.lucro.de/en/>





Philip

Art. no.	Width ex stock	Delivery time four weeks	Size
21180	M, W	S, +10, Upper	37-50



Björn

Art. no.	Width ex stock	Delivery time four weeks	Size
21030	M, W	S, +10, Upper	37-50



Jens

Art. no.	Width ex stock	Delivery time four weeks	Size
21260	M, W, +10	S, Upper	37-50



Heinrich

Art. no.	Width ex stock	Delivery time four weeks	Size
21399	S, M, W, +10	Upper	40-47



# kinetic



## The slimmer one

The LucRo kinetic is the slimmer-cut alternative that also provides ulcer-recurrence protection.

The feature of the LucRo kinetic line is a modified fit. We have given the last form a slimmer and narrower rear foot in the heel. This gives the shoe a different look and ensures an optimised heel guide especially for slimmer heels.

### Perfect all-round protection plus:

- Last in a narrower natural form
- Special kineticBags in the soles that relieve the joints



### Brigitte

Art. no.	Width ex stock	Delivery time four weeks	Size
25020	M, W	S, +10, Upper	36-43

## Women



### Elke

Art. no.	Width ex stock	Delivery time four weeks	Size
25112	M, W	S, +10, Upper	36-43



### Sabine

Art. no.	Width ex stock	Delivery time four weeks	Size
25070	S, M, W, +10	Upper	34-45



### Sabine

Art. no.	Width ex stock	Delivery time four weeks	Size
25375	M	S, W, +10, Upper	36-43



### Angelika

Art. no.	Width ex stock	Delivery time four weeks	Size
25413	M, W, +10	S, Upper	36-43

The models shown in this brochure are a choice of our collection.  
In order to see Schein's complete shoe range please go to: <https://www.lucro.de/en/>



Klaus

Art. no.	Width ex stock	Delivery time four weeks	Size
24020	M, W	S, +10, Upper	40–47



Rolf

Art. no.	Width ex stock	Delivery time four weeks	Size
24052	M, W, +10	S, Upper	37–50



Rolf

Art. no.	Width ex stock	Delivery time four weeks	Size
24053	M, W, +10	S, Upper	40–47



Rolf

Art. no.	Width ex stock	Delivery time four weeks	Size
24050	S, M, W, +10	Upper	37–50



Holger

Art. no.	Width ex stock	Delivery time four weeks	Size
24030	S, M, W, +10	Upper	37–50



Holger

Art. no.	Width ex stock	Delivery time four weeks	Size
24033	M	S, W, +10, Upper	40–47



Norbert

Art. no.	Width ex stock	Delivery time four weeks	Size
24040	M, W	S, +10, Upper	40–47



Norbert

Art. no.	Width ex stock	Delivery time four weeks	Size
24043	M, W	S, +10, Upper	40–47

The models shown in this brochure are a choice of our collection.  
 In order to see Schein's complete shoe range please go to: <https://www.lucro.de/en/>

## The sporty one

**New to the range: the specialist shoe for people with diabetes who want more active leisure time.**

The latest LucRo collection also has the protective properties for diabetes-friendly foot care.

### Women



Anke

Art. no.	Width ex stock	Delivery time four weeks	Size
30010	S, M, W, +10	Upper	36-43

Thanks to its anti-slip sole and sturdy design, it ensures stability especially over rough terrain.

### Perfect all-round protection plus:

- Last in a narrower natural form
- Anti-slip, sturdy TPU outsole
- Innovative, joint-protecting midsole



Anke

Art. no.	Width ex stock	Delivery time four weeks	Size
30015	M, W	S, +10, Upper	36-43



Includes spare pair of laces: LucRo flat laces in blue, 130 cm

Maren

Art. no.	Width ex stock	Delivery time four weeks	Size
30028	S, M, W, +10	Upper	36-43



Maren

Art. no.	Width ex stock	Delivery time four weeks	Size
30023	M, W	S, +10, Upper	36-43



Nina

Art. no.	Width ex stock	Delivery time four weeks	Size
30030	S, M, W, +10	Upper	36-43



Bianka

Art. no.	Width ex stock	Delivery time four weeks	Size
30062	M, W, +10	S, Upper	36-43

The models shown in this brochure are a choice of our collection. In order to see Schein's complete shoe range please go to: <https://www.lucro.de/en/>



Mario

Art. no.	Width ex stock	Delivery time four weeks	Size
28000	S, M, W, +10	Upper	40-50



Mario

Art. no.	Width ex stock	Delivery time four weeks	Size
28009	M, W, +10	S, Upper	40-50



Peco

Art. no.	Width ex stock	Delivery time four weeks	Size
28012	M	S, W, +10, Upper	40-47



Includes spare pair of laces: LucRo flat laces in dark blue, 140 cm

Peco

Art. no.	Width ex stock	Delivery time four weeks	Size
28018	M, W	S, +10, Upper	40-47



Manuel

Art. no.	Width ex stock	Delivery time four weeks	Size
28020	M, W	S, +10, Upper	40-47



Manuel

Art. no.	Width ex stock	Delivery time four weeks	Size
28024	S, M, W, +10	Upper	40-50



Jochen

Art. no.	Width ex stock	Delivery time four weeks	Size
28042	M, W, +10	S, Upper	40-50



Jochen

Art. no.	Width ex stock	Delivery time four weeks	Size
28048	M, W	S, +10, Upper	40-47

The models shown in this brochure are a choice of our collection.  
In order to see Schein's complete shoe range please go to: <https://www.lucro.de/en/>

## The adaptable one

In addition to proven ulcer-recurrence protection, the LucRo ergonic also provides ideal features for rheumatism care.

A recent diabetes study shows that the LucRo ergonic effectively protects against recurring foot ulcers. Its joint-friendly features also reduce the symptoms of rheumatism patients.

This is why the LucRo specialist shoe is listed for both indications.

### Perfect all-round protection plus:

- Last in a natural form
- Stiffening sole for optimum pressure distribution (must be firmly glued in for diabetes!)

## Women



Gisela

Art. no.	Width ex stock	Delivery time four weeks	Size
27010	M, W	S, Upper	36–43



Gisela

Art. no.	Width ex stock	Delivery time four weeks	Size
27011	M	S, W, +10, Upper	36–43



Bärbel

Art. no.	Width ex stock	Delivery time four weeks	Size
27049	M, +10	S, W, Upper	36–43



Bärbel

Art. no.	Width ex stock	Delivery time four weeks	Size
27042	S, M, W, +10	Upper	34–45



Marianne

Art. no.	Width ex stock	Delivery time four weeks	Size
27070	M	S, W, +10	36–43



Carola

Art. no.	Width ex stock	Delivery time four weeks	Size
27050	M, W, +10	S, Upper	36–43

The models shown in this brochure are a choice of our collection. In order to see Schein's complete shoe range please go to: <https://www.lucro.de/en/>



Lothar

Art. no.	Width ex stock	Delivery time four weeks	Size
26019	M, W	S, Upper	40-47



Lothar

Art. no.	Width ex stock	Delivery time four weeks	Size
26010	M, W	S, Upper	40-47



Jürgen

Art. no.	Width ex stock	Delivery time four weeks	Size
26020	M, W, +10	S, Upper	37-50



Jürgen

Art. no.	Width ex stock	Delivery time four weeks	Size
26023	M, W, +10	S, Upper	37-50



Karl

Art. no.	Width ex stock	Delivery time four weeks	Size
26043	M, W	S, +10, Upper	40-47



Reiner

Art. no.	Width ex stock	Delivery time four weeks	Size
26070	M, W	S, +10	37-50



Manfred

Art. no.	Width ex stock	Delivery time four weeks	Size
26032	M	S, W, +10, Upper	40-47



Detlef

Art. no.	Width ex stock	Delivery time four weeks	Size
26060	M, W, +10	S, Upper	40-47




The models shown in this brochure are a choice of our collection.  
In order to see Schein's complete shoe range please go to: <https://www.lucro.de/en/>



**Active and flexible –**  
without any strain or pain.

In comfort shoes the most important criteria are that the footwear must be not only comfortable but also functional. This is why our feel-good shoes are available in various sizes, widths and with varying heel heights. In this way, we can meet your individual needs for foot-friendly, comfortable shoes.

The technology behind actiflex: the stretch leather in the forefoot area is ideal for foot problems such as hallux valgus and hammer toe. It is stretchy and padded in such a way that the shoe does not look wide, but exudes natural elegance—for a beautiful and relaxed feel while walking.



**Our actiflex comfort shoes are particularly suitable for:**

- People who have to wear insoles
- Patients with mild to moderate foot conditions such as bunions, hammer toes and claw toes, as well as anyone who appreciates good shoes
- All people who set great store by pressure-free wearing comfort

**Our actiflex comfort shoes are characterised by:**

- Fashionable designs
- The greatest possible flexibility
- The use of first-class materials
- Continuous quality controls

**Explanations of symbols**



**Hallux Valgus**



**Hammer and claw toes**



**'all-in-one'**  
Hallux Valgus /  
Hammer and claw toes







Art. no.	Width	suitable for	Size
087016	G	Hallux Valgus	3-9



Art. no.	Width	suitable for	Size
084080	G	Hallux Valgus	3-9



Art. no.	Width	suitable for	Size
083101	G	'all-in-one'	2½-9



Art. no.	Width	suitable for	Size
080011	H	'all-in-one'	2½-9



Art. no.	Width	suitable for	Size
083141	H	'all-in-one'	2½-9



Art. no.	Width	suitable for	Size
087073	H	Hallux Valgus	3-9



Art. no.	Width	suitable for	Size
083101	H	Hallux Valgus	2½-9



Art. no.	Width	suitable for	Size
089095	H	Hammer and claw toes	3-9

The models shown in this brochure are a choice of our collection.  
 In order to see Schein's complete shoe range please go to: <https://actiflex.de/en/home.html>

## Women



Art. no.	Width	suitable for	Size
089096	H	Hammer and claw toes	3-9



Art. no.	Width	suitable for	Size
081552	K	Hammer and claw toes	3-9



Art. no.	Width	suitable for	Size
089509	K	Hammer and claw toes	3-9



Art. no.	Width	suitable for	Size
089523	K	Hallux Valgus	3-9



Art. no.	Width	suitable for	Size
087309	H	Hallux Valgus	3-9

The models shown in this brochure are a choice of our collection.  
 In order to see Schein's complete shoe range please go to: <https://actiflex.de/en/home.html>



Art. no.	Width	suitable for	Size
088003	H	Hammer and claw toes	6-12

Art. no.	Width	suitable for	Size
088710	H	Hammer and claw toes	6-12



Art. no.	Width	suitable for	Size
088009	H	'all-in-one'	6-12




## Comfortable and stylish

Supreme comfort doesn't come at the expense of fashionable design. This is proven by the models in our Ortho Lady & Gent collection, which meet the expectations of a modern, health-conscious generation.

In the ranges we have a wide selection of comfortable shoes which support the natural physical motion and counteract tiredness of the feet. It goes without saying that all shoes are suitable for use with loose insoles. With lovingly created details and refined appliquéés we have created a varied look combined with a high degree of comfort.

## Women

Around the inner ball area of some models we have incorporated a special material which provides relief in the area of the hallux valgus in case of deformations. These models are marked with this symbol .


### Comfort footbed

- Support the natural gait
- Excellent fit: good comfort
- Shock-absorbing
- Eudermic materials
- Breathable leather cover
- Low maintenance



### The structure

- 1 leather
- 2 EVA foam
- 3 Shock Absorber
- 4 hook and loop fastener



The models shown in this brochure are a choice of our collection.

In order to see Schein's complete shoe range please go to: <https://schein.de/en/produkte/schuhe.html?marke=olog>



Art. no.	Size
385002	2 1/2-8 1/2



Art. no.	Size
381080	2 1/2-8 1/2



Art. no.	Size
381371	2 1/2-8 1/2



Art. no.	Size
388291	2 1/2-8 1/2



Art. no.	Size
387003	2 1/2-8 1/2



Art. no.	Size
385266	2 1/2-8 1/2




Art. no.	Size
389820	2 1/2-8 1/2



Art. no.	Size
381280	2 1/2-8 1/2



## Women



Art. no.  
381051



Size  
2 1/2-8 1/2



Art. no.  
389700



Size  
2 1/2-8 1/2



Art. no.  
380504



Size  
2 1/2-8 1/2



Art. no.  
385553



Size  
2 1/2-8 1/2





Art. no.



Size

381653

2 1/2-8 1/2



Art. no.



Size

381683

2 1/2-8 1/2



Art. no.

Size

399000

6-12



Art. no.

Size

394013

6-12

Men



### The structure

- ① leather
- ② natural cork

# Post op shoes

## Perfect care— after surgery



The foot needs special protection after surgery. No pressure should be applied to the wound and the shoe should not rub on the foot at any point. The foot should also be stabilised and immobilised.

The post-op shoes shown offer different characteristics for the special requirements of the operated foot.

The following models are available:

- Post op shoe
- Forefoot-relief shoe
- Heel-relief shoe

Shoe size	Internal shoe length a	Internal shoe width b
33/34	224 mm	82 mm
35/36	236 mm	86 mm
37/38	250 mm	90 mm
39/40	264 mm	94 mm
41/42	278 mm	98 mm
43/44	292 mm	102 mm
45/46	306 mm	106 mm
47/48	320 mm	110 mm



### Indications

- Post-operative foot surgery
- Post-traumatic forefoot injuries
- Open wounds and ulcers
- Feet with strong tendencies to swelling



Art. no.	Model	Size
322146	Post op shoe (short-term use), closed heel	33–48





### Indications



- Post-operative foot surgery
- Post-traumatic forefoot injuries
- Open wounds and ulcers
- Feet with strong tendencies to swelling



Art. no.	Model	Size
322046	Post op shoe (short-term use), closed heel	33/34–47/48



### Indications



- Post-operative forefoot surgery
- Post-traumatic forefoot injuries
- Open wounds and ulcers located at the dorsal foot or toes
- Acute swelling of the forefoot due to bursitis or gout-attacks



Art. no.	Model	Size
322147	Forefoot-relief shoe	33/34–47/48



### Indications



- Post-operative forefoot surgery
- Post-traumatic forefoot injuries
- Open wounds and ulcers located at the dorsal foot or toes
- Acute swelling of the forefoot due to bursitis or gout-attacks



Art. no.	Model	Size
322047	Forefoot-relief shoe	33/34–47/48

# Post op shoes



## Indications



- Post-operative foot surgery
- Post-traumatic forefoot and back foot injuries
- Open wounds and ulcers
- Feet with strong tendencies to swelling



Art. no.	Model	Size
322048	Post op shoe (short-term use), open heel	33/34–47/48



## Indications



- Post-operative hindfoot surgery
- Post-traumatic hindfoot injuries
- Open wounds and ulcers in the heel area
- With painful and infective illnesses of the heel area



Art. no.	Model	Size
322003	Heel-relief shoe	33/34–47/48

## Suitable for:

- 322146
- 322147




## Suitable for:

- 322046
- 322047
- 322048
- 322003



Art. no.	Model
322149	Vamp for post op shoes

Art. no.	Model
322049	Vamp for Post-Op shoes



# Novaped

by schein

## The basis for improved performance

People who are very active in their free time often feel significant strain on their feet.

The aim of NovaPED performance is to counteract stress and discomfort day to day or during sporting activities.

Three levels, one aim:

the **multi-zone level** targets the ressure zones of the foot.



The **functional level** is made up of three materials, each with a different function:



**SOFT**

Soft cushioning

Multiform hybrid



**MEDIUM**

Shock-absorbing

p<sup>2</sup> absorb



**PLUS**

Springy, energy-returning

p<sup>2</sup> rebound



The **insole level**, which is meiculously pre-shaped, has a slightly pronating effect on step initiation as a result of the proportion of approx. 30° (grey) and approx. 40° (black) Shore A hardness.

Optionally available with a perforated, breathable OnDry cover.

The blank can be reshaped at 100°C.



Art. no.	Description	Size
978110000	NovaPED performance SOFT	35–48
978111000	NovaPED performance SOFT, without metatarsal pad	35–48



Art. no.	Description	Size
978110001	NovaPED performance SOFT with cover OnDry black perforated	35–48
978111001	NovaPED performance SOFT with cover OnDry black perforated without metatarsal pad	35–48



Art. no.	Description	Size
978120000	NovaPED performance MEDIUM	35–48
978121000	NovaPED performance MEDIUM, without metatarsal pad	35–48



Art. no.	Description	Size
978120001	NovaPED performance MEDIUM with cover OnDry black perforated	35–48
978121001	NovaPED performance MEDIUM with cover OnDry black perforated without metatarsal pad	35–48



Art. no.	Description	Size
978130000	NovaPED performance PLUS	35–48
978131000	NovaPED performance PLUS, without metatarsal pad	35–48



Art. no.	Description	Size
978130001	NovaPED performance PLUS with cover OnDry black perforated	35–48
978131001	NovaPED performance PLUS with cover OnDry black perforated without metatarsal pad	35–48



The orthopaedic insole blanks shown in this brochure are a choice of your range.  
In order to see Schein's complete program please go to: <https://www.schein.de/en/download/prospekte.html>

## CARBON EDITION

The high level of functionality of the S90 Carbon insole is down to the lightness, stability and low structural height of carbon as a material. The material particularly demonstrates its strengths in treatment of plantar fasciitis and hallux rigidus. The S90 film protects the carbon core and allows the customer to see the sophisticated design.

Novaped S90 Carbon insoles are pre-shaped and available in a glued version and as a base (chassis). The chassis can be built up with the materials of your choice and thus individually adjusted to different foot types.

### Hallux Rigidus



Without cover



1.5 mm padding

Art. no. left	Art. no. right	Padding	Cover	Size
983220000	983120000	–	–	35–48
983221000	983121000	p <sup>2</sup> , ca. 1.5 mm padding, toe area thinner	Microfibre OnDry by onSteam, black	35–48

### Compensation insole



Without cover



1.5-mm padding

Art. no. left	Art. no. right	Padding	Cover	Size
983240000	983140000	–	–	35–48
983241000	983141000	Continuous p <sup>2</sup> padding, approx. 1.5 mm	Microfibre OnDry by onSteam, black	35–48





## Carbon-soles

### Items in stock

Art. no. pair	Description flat	Size
0479000	ladies, thickness: 1,6 mm	36–41
0479001	men, thickness: 1,6 mm	41–47
0479050	ladies, thickness: 1,9 mm	36–41
0479051	men, thickness: 1,9 mm	41–47

Art. no. pair	Description preformed	Size
954795000	ladies, thickness: 1,6 mm	36–41
954796000	men, thickness: 1,6 mm	41–47
954797000	ladies, thickness: 1,9 mm	36–41
954798000	men, thickness: 1,9 mm	41–47

## Plantar fasciitis

### Plantar fasciitis insole with heel spur cut-out

Without  
cover



With cover/  
padded cover/  
foot-side relief pad



Art. no. pair	Padding	Cover	Size
982810000	–	–	35–48
982811000	p <sup>2</sup> -relief pad 2 mm, continuous p <sup>2</sup> padding 1.5 mm	Microfibre OnDry by onSteam, black	35–48

### Plantar fasciitis insole with heel spur cut-out, tendon cut-out and padded ring on shoe side

Without cover/  
heel pad



With cover/  
heel pad/  
padded cover



Art. no. pair	Padding	Cover	Size
982800000	Heel pad	–	35–48
982801000	Heel pad, continuous p <sup>2</sup> padding 2 mm	Microfibre OnDry by onSteam, black	35–48

The orthopaedic insole blanks shown in this brochure are a choice of your range.

In order to see Schein's complete program please go to: <https://www.schein.de/en/download/prospekte.html>

## sports

Running / Walking /  
Nordic Walking



3.5 mm padding

Art. no. pair	Padding	Cover	Size
98250000	Continuous p <sup>2</sup> /EVA sandwich padding, approx. 3.5 mm	Microfibre OnDry by onSteam, black	35–48

## Football



3.5 mm padding

Art. no. pair	Padding	Cover	Size
982501000	Continuous p <sup>2</sup> /EVA sandwich padding, approx. 3.5 mm	Microfibre OnDry by onSteam, black	35–48

## Cycling



2.0 mm padding

Art. no. pair	Padding	Cover	Size
982502000	Continuous p <sup>2</sup> padding, approx. 2 mm	Microfibre OnDry by onSteam, black	35–48




### Corrective/supportive

We offer shaping elements lined with a heat-activatable adhesive film for individual adjustment of our Novaped S90 insoles. These elements can be attached quickly, safely and permanently.

Thanks to its high Shore A hardness of approx. 58°, ARU®leicht is particularly suitable for corrective use.

As well as the sheet of material, we also offer you prefabricated elements that you can use to quickly treat various indications.



Art. no.	Description
964500864	Sheet, approx. 190 x 390 x 4 mm
964500866	Sheet, approx. 190 x 390 x 6 mm



Three usage possibilities with just one element.

Art. no.	Description		Size
964520866	Heel wedge 6 mm	1	35–39
		2	40–43
		3	44–48

Art. no.	Description		Size
964510864	Pronation/Supination 4 mm	1	35–39
		2	40–43
		3	44–48

The orthopaedic insole blanks shown in this brochure are a choice of your range. In order to see Schein's complete program please go to: <https://www.schein.de/en/download/prospekte.html>

## Orthopaedic foot supports with heel plateau

### Effective treatment

The Novaped plateau orthopaedic foot supports offer maximum stability together with optimal comfort.


With its extended support surface, the heel plateau provides stability and improves transmission of force onto the sole of the shoe. The heel plateau also features an integrated rolling edge that supports controlled step initiation, thus reducing the speed of pronation.

Novaped S90 cupped insoles with heel plateau stand out by stabilising the rearfoot thanks to an enlarged support surface.

Pronation and supination corrections can be made directly on the insole. Specific points on the S90 insoles can be reshaped very effectively with a hot-air gun.

### Benefits


- Round on the inside, straight on the outside
  - Appropriate for the anatomy of the heel on the inside, with an extended support surface on the outside
- Freedom to adjust
  - For pronation and supination corrections
- Controlled step initiation
  - Thanks to the tapered heel impact zone
- Stabilised tread
  - Thanks to extended support surface on the heel base



With appropriate further adjustment by a healthcare technician, the heel plateau can stabilise heel impact and optimise gait.


## Round on the inside

The rounded heel and the cupped design provide a perfect fit for the foot. This is beneficial for pressure distribution and ensures the desired treatment result while also being comfortable to wear.



## Straight on the outside

With its extended support surface, the heel plateau stops the insole from rocking and, at the same time, improves transmission of force from the insole onto the sole of the shoe.




## Freedom to adjust

Pronation and supination corrections can be made directly on the plateau (without the need for additional material).

## Controlled step initiation

The integrated rolling edge aids controlled step initiation.




## Supporting orthopaedic insoles with heel plateau

They are manufactured using the tried-and-tested 'film technique'. The special sandwich construction gives the insoles a low mounting height, flexible stability and soft, stable edges.

The thin construction means the finished plateau insoles will fit into ready-made shoes with ease.

Thermal mouldability at 100°C allows additional customisation by healthcare technicians.



Art. no.	Description	Size
959310000	Novaped plateau base	35–48
959312000	Novaped plateau, 2 mm p <sup>2</sup> resorb padding and perforated OnDry cover	35–48

The orthopaedic insole blanks shown in this brochure are examples of our range. In order to see Schein's complete range, please go to: <https://www.schein.de/en/download/prospekte.html>

## Cushioning orthopaedic foot supports with heel plateau

They are manufactured with a multi-zone technique on the insole base. The rearfoot is made of a cork composite and provides stability, while the shock-absorbing forefoot

Multiforte provides relief for the metatarsals. The heel plateau integrated into the rearfoot enables pronation and supination corrections.

### Padding

- Multiform
- Thickness 2 mm
- Beige
- Approx. 30° Shore A
- For shock absorption on the plantar surface of the foot

### Cover

- Perforated OnDry cover to stabilise, reinforced to 1.1 mm
- Beige
- Moisture-regulating



### Heel plateau

- For improved adjustment to the insole
- Allows targeted further adjustment to support pronation or supination



- Shock-absorbing forefoot component
- Multiform
- Grey
- Approx. 30° Shore A
- Stabilising rear-foot component
- Multikork
- Beige
- Approx. 60° Shore A

### Padding

- Multiform
- Thickness 2 mm
- Grey
- Approx. 30° Shore A
- For shock absorption on the plantar surface of the foot

### Cover

- Perforated OnDry cover
- Black
- Moisture-regulating



### Heel plateau

- For improved adjustment to the insole
- Allows targeted further adjustment to support pronation or supination



- Shock-absorbing forefoot component
- Multiform
- Grey
- Approx. 30° Shore A
- Stabilising rearfoot component
- Multiforte
- Black
- Approx. 40° Shore A


Art. no.	Description	Size
959322000	Novaped plateau with divided base Multikork/Multiform, 2 mm beige Multiform, reinforced perforated OnDry cover	35–48
979332000	Novaped plateau with divided base (Multiform/Multiforte), 2 mm grey Multiform, perforated OnDry cover	35–48

## Cupped orthopaedic foot supports with heel plateau

They offer a high level of comfort and support, and can correct foot malpositions with professional further adjustment.

For optimal treatment results, it is important that the foot, insole and shoe are perfectly harmonised.

The orthopaedic foot support can be thermally moulded at 100°C.



The orthopaedic insole blanks shown in this brochure are examples of our range.

In order to see Schein's complete range, please go to: <https://www.schein.de/en/download/prospekte.html>

- **Medium** version
- Without metatarsal pad

Art. no.	Description	Size
980189000	Novaped plateau with PiGreco cover, cupped shape	19/20–48



- **Narrow** version
- Slightly pressed metatarsal pad

Art. no.	Description	Size
980360000	Novaped plateau with 2 mm Multiform and PiGreco cover, cupped shape	19/20–47



- **Wide** version
- Slight retrocapital support

Art. no.	Description	Size
980260000	Novaped plateau with 2 mm Multiform, cupped shape	19/20–42





# Novatherm insoles

Novatherm forming sets or insoles are highly suitable for adjustment to the needs of the customer.

The forming set, made of thermoplastic materials, is best shaped with a vacuum press at a temperature of 80–90°C. (We recommend setting the temperature of the oven or heating plate to 150°C.)



The finished insole blank can be heated and reshaped at specific points with the help of a hot-air gun.

The cooled blank is then very suitable for grinding or other further adjustment and is characterised by a high level of stability.

## Novatherm forming sets

Art. no.	Description	Size
986089000	Microfibre children's set, top cover	19/20–33/34, 3–48

Art. no.	Description	Size
986020000	Adult set, small	35–42

Art. no.	Description	Size
986030000	Set, wide	19/20–48



## Novatherm insoles

- Full length and width

Art. no.	Description	Size
985146454	'Comfort long'	27/28–48



- With especially high and stable edges

Art. no.	Description	Size
981289000	'Lange'	19/20–48



- Full length and width
- With especially high and stable edges

Art. no.	Description	Size
986543454	'Aktiv'	35–48



- Full length and width
- With especially high and stable edges
- Top layer PiGreco, black
- Thin, long PU 'Aktiv with PU' padding (1.5 mm)

Art. no.	Description	Size
986539000	'Aktiv with PU'	35–48




The orthopaedic insole blanks shown in this brochure are examples of our range. In order to see Schein's complete range, please go to: <https://www.schein.de/en/download/prospekte.html>

# Novaped moulding sets

## Customised insoles

- Suitable functional materials are selected for the supportive base in accordance with the type of treatment required. Materials with cushioning, padding or shock-absorbing effects are available.
- The OnDry cover with hydrophilic and antibacterial properties protects the materials below and ensures the foot is comfortable.
- All materials are carefully aligned and harmonised, meaning that no additional adhesive is required.
- The insoles have to be individually produced in a vacuum process.



- Narrow/wide base

Art. no.	Description	Size
586040000	Narrow base	35-52
586041000	Wide base	35-52




- Padding materials

Art. no.	Description	Size
464705700	p <sup>2</sup> resorb/green, wide, 2 mm	35-52
464705701	p <sup>2</sup> resorb/red, wide, 2 mm	35-52
464705702	p <sup>2</sup> resorb/yellow, wide, 2 mm	35-52
464705712	p <sup>2</sup> resorb/black, wide, 2 mm	35-52



- Cover materials

Art. no.	Description	Size
464705542	OnDry cover/gobi dark, wide	35-52
464705543	OnDry cover/grey, wide	35-52
464705544	OnDry cover/black, wide	35-52



The orthopaedic insole blanks shown in this brochure are examples of our range.  
In order to see Schein's complete range, please go to: <https://www.schein.de/en/download/prospekte.html>

The Novaped sports foot support blanks are designed for different sports with their specific movement and strain patterns. They ensure high wearing comfort, prevent excessive strain and contribute to the alleviation of complaints. The blanks are pre-shaped.

The Novaped sports blanks, made of thermo-plastic materials, are best shaped with a vacuum press at a temperature of 80–90°C. (We recommend setting the temperature of the oven or heating plate to 150°C.)

Alternatively, the pre-shaped blank can be heated and reshaped at specific points with the help of a hot-air gun.



The orthopaedic insole blanks shown in this brochure are examples of our range. In order to see Schein's complete range, please go to: <https://www.schein.de/en/download/prospekte.html>

Art. no.	Description	Size
987620000	full length, wide, natural metatarsal pad	35-53
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	Microfibre, orange

Running



Art. no.	Description	Size
987624000	full length, narrow, natural metatarsal pad	35-46
Base	Top layer	Top Size
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	Microfibre, orange

Running



Art. no.	Description	Size
987627000	full length, wide, natural metatarsal pad, raised outer edge	35-53
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	Microfibre, orange

Running

Special support with raised under edge



Art. no.	Description	Size
987626000	full length, narrow, natural metatarsal pad, raised outer edge	35-46
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	Microfibre, orange

Running

Special support with raised under edge



Running  
Special support with cut-out of the plantar fascia tendon



Art. no.	Description	Size
987623000	full length, wide, natural metatarsal pad, cut-out of the plantar fascia	35–49
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	Microfibre, orange

Running  
Special support  
Construction kit



Art. no.	Description	Size
987621000	full length, wide, natural metatarsal pad, cut-out of the plantar fascia	35–49
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm, loose pad made of multiform, white, Shore A approx. 30°	

Football



Art. no.	Description	Size
987610000	full length, medium width, natural metatarsal pad from size 35	29/30–48
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	ARU X-Static, comb pattern

Ball sports



Art. no.	Description	Size
987661000	full length, wide, natural metatarsal pad	35–53
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	Microfibre, blue



Art. no.	Description	Size
987668000	full length, narrow, natural metatarsal pad	35–46
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	Microfibre, blue

## Ball sports



Art. no.	Description	Size
987616000	full length, wide, natural metatarsal pad	35–53
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	Microfibre, green

## Racket



Art. no.	Description	Size
987615000	full length, narrow, natural metatarsal pad	35–46
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	Microfibre, green

## Racket



Art. no.	Description	Size
987686000	full length, medium width	35–48
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 1,5 mm	Microfibre, grey

## Bike



The orthopaedic insole blanks shown in this brochure are examples of our range.

In order to see Schein's complete range, please go to: <https://www.schein.de/en/download/prospekte.html>



## Golf



Art. no.	Description	Size
987618000	full length, wide, natural metatarsal pad	35–53
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	Microfibre, beige

## Golf



Art. no.	Description	Size
987617000	full length, narrow, natural metatarsal pad	35–46
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	Microfibre, beige

## Skating



Art. no.	Description	Size
987632000	full length, wide	35–48
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	ARU X-Static, comb pattern

## Skating



Art. no.	Description	Size
987631000	full length, narrow	35–46
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	ARU X-Static, comb pattern

Art. no.	Description	Size
987657000	full length, wide	35–48
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	ARU X-Static, comb pattern

## Cross-country skiing



Art. no.	Description	Size
987656000	full length, narrow	35–46
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	ARU X-Static, comb pattern

## Cross-country skiing



Art. no.	Description	Size
987635000	full length, medium width, natural metatarsal pad	35– 48
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 3 mm	ARU X-Static, comb pattern

## Downhill skiing



Art. no.	Description	Size
987638000	full length, medium width, natural metatarsal pad	35–48
Base	Top layer	Top cover
Sandwich construction, transparent, PE reinforcement	p <sup>2</sup> , skin colour, Shore A approx. 18°, 2 mm	ARU X-Static, comb pattern

## Snowboard



The orthopaedic insole blanks shown in this brochure are examples of our range.

In order to see Schein's complete range, please go to: <https://www.schein.de/en/download/prospekte.html>

Novaped soft insoles are primarily made of Multiform (approx. 30° Shore A) and Multiforte (approx. 40° Shore A) materials, rounded off with stiffening (THK) or cushioning (PPT or similar) elements.

Some models are optionally available with OnDry (microfibre) covers.

Novaped soft insoles are very suitable for grinding, thermal reshaping and adhering to other materials.



Art. no.	Description	Size
972030300	Short/ medium width	35–48
<b>Base</b>	<b>Top layer</b>	<b>Top cover</b>
Multiform, white, smooth, Shore A approx. 30°, 15 mm		



Art. no.	Description	Size
971320000	Full length, extra wide, natural metatarsal pad	35–48
<b>Base</b>	<b>Top layer</b>	<b>Top cover</b>
Multiform, white, Shore A approx. 30°, 15 mm		



Art. no.	Description	Size
977013300	Quarter cup, full length, no meta-pad imprint	19/20–48
<b>Base</b>	<b>Top layer</b>	<b>Top cover</b>
Multiform, white, Shore A approx. 30°, 10 mm: size 19/20–33/34 15 mm: size 35–48		

Art. no.	Description	Size
977043000	Quarter cup, full length, natural metatarsal pad (adult sizes only)	19/20–48
Base	Top layer	Top cover
Multiforte, white, smooth, Shore A approx. 40°, 12 mm: size 19/20–33/34 15 mm: size 35–48		



Art. no.	Description	Size
971306300	Half cup, full length, no meta-pad imprint	35–48
Base	Top layer	Top cover
Multiform, white, smooth, Shore A approx. 30°, 12 mm: size 19/20–33/34 15 mm: size 35–48		



Art. no.	Description	Size
477014000	Half shell, full length, natural metatarsal pad	35–47
Base	Top layer	Top cover
Multiform, white, smooth, Shore A approx. 30°, 10 mm: size 19/20–33/34 15 mm: size 35–47		



The orthopaedic insole blanks shown in this brochure are examples of our range.  
In order to see Schein's complete range, please go to: <https://www.schein.de/en/download/prospekte.html>



Art. no.	Description	Size
470610000	Full length, medium width, natural metatarsal pad (adult sizes only)	19/20–48
<b>Base</b>	<b>Top layer</b>	<b>Top cover</b>
Multiforme, white, Shore A approx. 40°, 10 mm, THK reinforcement	Multiform, white, Shore A approx. 30°, 4 mm, perforated	



Art. no.	Description	Size
977753000	Natural shape, full length, natural metatarsal pad	35–48
<b>Base</b>	<b>Top layer</b>	<b>Top cover</b>
Multiform, white, Shore A approx. 30°, 10 mm, THK reinforcement	Multiform, white, Shore A approx. 30°, 7 mm	



Art. no.	Description	Size
470143000	Half cup, full length, natural metatarsal pad	35–48
<b>Base</b>	<b>Top layer</b>	<b>Top cover</b>
Multiforme, white, Shore A approx. 40°, 10 mm	Multiform, white, Shore A approx. 30°, 10 mm	



Art. no.	Description	Size
470141000	Half cup, full length, natural metatarsal pad (adult sizes only)	19/20–48
<b>Base</b>	<b>Top layer</b>	<b>Top cover</b>
Multiforme, white, Shore A approx. 40°, 10 mm	Multiforme, white, Shore A approx. 40°, 4 mm	

Art. no.	Description	Size
470427209	Full length, narrow, natural metatarsal pad (adult sizes only)	19/20–48

Base	Top layer	Top cover
Multiform, white, Shore A approx. 30°, 7 mm, THK reinforcement	Multiform, blue, Shore A approx. 30°, 4 mm, perforated	



Art. no.	Description	Size
973326000	Full length, wide, low cup, natural metatarsal pad (adult sizes only)	19/20–48

Base	Top layer	Top cover
Multiforte, white, Shore A approx. 40°, 10 mm, THK reinforcement	Multiform, blue, Shore A approx. 30°, 4 mm, perforated	



Art. no.	Description	Size
977713000	Full length, narrow, quarter cup, natural metatarsal pad (adult sizes only)	19/20–48

Base	Top layer	Top cover
Multiform, white, Shore A approx. 30°, 10 mm, THK reinforcement	Multiform, blue, Shore A approx. 30°, 4 mm, perforated	



Art. no.	Description	Size
470374000	half cup, full length, natural metatarsal pad	35–48

Base	Top layer	Top cover
Multiforte, white, Shore A approx. 40°, 20 mm	Multiform, blue, Shore A approx. 30°, 4 mm, perforated	



The orthopaedic insole blanks shown in this brochure are examples of our range.

In order to see Schein's complete range, please go to: <https://www.schein.de/en/download/prospekte.html>



Art. no.	Description	Size
470899000	Natural shape, full length, natural metatarsal pad (adult sizes only)	19/20–50
Base	Top layer	Top cover
Multiforte, blue, Shore A approx. 40°, 5 mm, LA-wedge, THK reinforcement	Multicom, blue/blue, Shore A approx. 30°/21°, 4 mm, perforated	



Art. no.	Description	Size
471646000	Full length, medium width	35–50
Base	Top layer	Top cover
Multiforte, blue, Shore A approx. 40°, 5 mm, THK reinforcement	Multicom, blue/blue, Shore A approx. 30°/21°, 4 mm, perforated	PiGreco by schein, black



Art. no.	Description	Size
474097000	Natural shape, full length, natural metatarsal pad	35–46
Base	Top layer	Top cover
Multiforte, blue, Shore A approx. 40°, 12 mm, LA-wedge, THK reinforcement	Multicom, blue/blue, Shore A approx. 30°/21°, 4 mm, perforated	



Art. no.	Description	Size
470807000	Cup shape, full length, natural metatarsal pad	19/20–47
Base	Top layer	Top cover
Multicork, Shore A approx. 60°, 10 mm, THK reinforcement	Multiform, blue, Shore A approx. 30°, 5 mm, perforated	



Art. no.	Description	Size
470432300	Full length, medium width	35–48

Base	Top layer	Top cover
Multiform, black, Shore A approx. 30°, 10 mm, THK reinforcement	Multiform, black, Shore A approx. 30°, 5 mm	



Art. no.	Description	Size
470432000	Full length, medium width, natural metatarsal pad	35–48

Base	Top layer	Top cover
Multiform, black, Shore A approx. 30°, 10 mm, THK reinforcement	Multiform, black, Shore A approx. 30°, 5 mm	



Art. no.	Description	Size
470434000	Full length, wide, natural metatarsal pad	35–48

Base	Top layer	Top cover
Multiform, black, Shore A approx. 30°, 10 mm, THK reinforcement	Multiform, black, Shore A approx. 30°, 5 mm	



Art. no.	Description	Size
470427000	Full length, narrow, natural metatarsal pad (adult sizes only)	19/20–48

Base	Top layer	Top cover
Multiform, white, Shore A approx. 30°, 7 mm, THK reinforcement	Multiform, pink, Shore A approx. 30°, 4 mm, perforated	




The orthopaedic insole blanks shown in this brochure are examples of our range.

In order to see Schein's complete range, please go to: <https://www.schein.de/en/download/prospekte.html>



Art. no.	Description	Size
470295000	Natural shape, full length, natural metatarsal pad	36–48
Base	Top layer	Top cover
Multiform, firm, white, Shore A approx. 50°, 15 mm	Lunairflex, skin colour, Shore A approx. 22°, 4 mm, perforated	



Art. no.	Description	Size
470770133	Full length, wide, natural metatarsal pad (adult sizes only)	19/20–53
Base	Top layer	Top cover
Multiform, white, Shore A approx. 30°, 10 mm		Neat's leather, perforated, soft



Art. no.	Description	Size
470600000	Cup shape, full length, natural metatarsal pad (adult sizes only)	19/20–48
Base	Top layer	Top cover
Thermaforte, blue, Shore A approx. 40°, 8 mm, THK reinforcement	Multiform, multicolor seven colours, Shore A approx. 30°, 4 mm	

The orthopaedic insole blanks shown in this brochure are examples of our range.  
 In order to see Schein's complete range, please go to: <https://www.schein.de/en/download/prospekte.html>





Our Novaped soft insoles for LucRo special shoes are multi-layer, pressure-re-distributing, soft-padded cushioning insoles. The insoles provide moderate longitudinal arch support and have no meta pad. The shape is adapted to our LucRo protective shoe models. The insoles are available with or without a cover.

#### Optional with

- Cover Scheinova
- Padding layer Lunairmed, Shore A approx. 18°, 6 mm
- Stiffening Flexible special core
- Foundation Lunairflex, Shore A approx. 22°, 6 mm

basic



Art. no.	Description		Size
972204...	LucRo classic	Women	34-45
972105...	LucRo classic	Men	37-50
972214...	LucRo kinetic LucRo sportic	Women	34-45
972115...	LucRo kinetic LucRo sportic	Men	37-50
979481...	LucRo ergonic	Women	34-45
979480...	LucRo ergonic	Men	37-50

-000 without cover

-637 with cover

comfort

Novaped soft is a multi-layer, pressure-re-distributing, soft-padded cushioning insole. The anatomical shape of the last consists of cupped heel guidance, a recess in the base V and retrocapital transverse support.

The insoles are available in women's and men's versions, either with or without a microfibre cover. The corresponding shapes are suitable for adaptation in our LucRo shoes for people with diabetes or rheumatism.

#### Optional with

Cover	Microfibre
Padding layer	Combination of: ARUfoam, Shore A approx. 14°, 3 mm p <sup>2</sup> resorb, Shore A approx. 10°, 2 mm ARUfant supersoft, Shore A approx. 25°, 4 mm
Stiffening	THK stiffening
Foundation	ARUfant, Shore A approx. 55°, 8 mm



Art. no.	Description	Size
971880000	Without cover Women	35–48
971881000	With cover Women	35–48
971870000	Without cover Men	35–52
971871000	With cover Men	35–52

The orthopaedic insole blanks shown in this brochure are examples of our range. In order to see Schein's complete range, please go to: <https://www.schein.de/en/download/prospekte.html>

# NovaPED blanks

Ideal for diabetic care

ARU® combi **TricoDirect**® materials are suitable for the production of individual footbeds. The combinations of three vulcanised material layers can be easily formed and further processed.

ARU® combi **TricoDirect**® stands out with high quality and good durability.

**EVA material** with three vulcanised layers of different Shore A.

ARU®**liso**, light cell EVA with very good resilience

ARU®**Reboflex**, EVA, closed-cell, high resetting features with GEL character.

Weight-specific material combinations, washable and disinfectable.



Heating recommendation (temperature/duration)



Cooling time recommendation (can be shortened with the application of a wet cloth)

## ARU®combi **TricoDirect**® 12


5 mm ARU® <b>liso</b> ,	approx. 20° Shore A
5 mm ARU® <b>Reboflex</b> ,	approx. 25° Shore A
12 mm ARU® <b>liso</b> ,	approx. 35° Shore A

## ARU®combi **TricoDirect**® 8


5 mm ARU® <b>liso</b> ,	approx. 20° Shore A
5 mm ARU® <b>liso</b> ,	approx. 35° Shore A
10 mm ARU® <b>liso</b> ,	approx. 45° Shore A

## ARU®combi **TricoDirect**® 10



5 mm ARU® <b>liso</b> ,	approx. 20° Shore A
5 mm ARU® <b>liso</b> ,	approx. 35° Shore A
10 mm ARU® <b>liso</b> ,	approx. 55/60° Shore A



Patients up to **70 kg**



Patients up to **90 kg**



Patients up to **105 kg**

Art. no.	Description
46 0005-212	Punching blank
103 403 12	Plate, 950 x 550 mm

Art. no.	Description
46 0005-208	Punching blank
103 403 08	Plate, 950 x 550 mm

Art. no.	Description
46 0005-210	Punching blank
103 403 10	Plate, 950 x 550 mm

The orthopaedic insole blanks shown in this brochure are examples of our range.

In order to see Schein's complete range, please go to: <https://www.schein.de/en/download/prospekte.html>







# Schein padding and bedding materials

## Multiform (sheet goods)

Polyethylene cellular material, multi-purposed cushioning material, ground both sides, closed-cell, extreme high resistance to pressure, very suitable for deep-drawing, kind to the skin, resistant to sweat, easy to process, washable,

thermoforming at: approx. 100°–130°C, hardness: approx. 30° Shore A

### application:

foot orthotics and footbeds, sport shoes



\* = perforated

Art. no.	Colour	Thickness approx. mm	Size approx. mm
010590010	white	1,0	1.100 x 1.100
010590020	white	2,0	1.100 x 1.100
010590030	white	3,0	1.100 x 1.100
010590040	white	4,0	1.100 x 1.100
010590050	white	5,0	1.100 x 1.100
010590060	white	6,0	1.100 x 1.100
010590070	white	7,0	1.100 x 1.100
010590100	white	10,0	1.100 x 1.100
010590150	white	15,0	1.100 x 1.100
010592020*	white	2,0	1.100 x 1.100
010592030*	white	3,0	1.100 x 1.100
010592040*	white	4,0	1.100 x 1.100
010592050*	white	5,0	1.100 x 1.100
010581070	green	7,0	1.100 x 1.100
010589020*	green	2,0	1.100 x 1.100
010589040*	green	4,0	1.100 x 1.100
010582020	blue	2,0	1.100 x 1.100
010582030	blue	3,0	1.100 x 1.100
010582040	blue	4,0	1.100 x 1.100
010582050	blue	5,0	1.100 x 1.100
010582060	blue	6,0	1.100 x 1.100
010582070	blue	7,0	1.100 x 1.100
010585020*	blue	2,0	1.100 x 1.100
010585030*	blue	3,0	1.100 x 1.100
010585040*	blue	4,0	1.100 x 1.100
010585050*	blue	5,0	1.100 x 1.100
010585060*	blue	6,0	1.100 x 1.100
010585070*	blue	7,0	1.100 x 1.100
010595020	skin coloured	2,0	1.100 x 1.100
010595030	skin coloured	3,0	1.100 x 1.100
010595040	skin coloured	4,0	1.100 x 1.100
010595050	skin coloured	5,0	1.100 x 1.100

Art. no.	Colour	Thickness approx. mm	Size approx. mm
010595060	skin coloured	6,0	1.100 x 1.100
010595070	skin coloured	7,0	1.100 x 1.100
010595080	skin coloured	8,0	1.100 x 1.100
010595100	skin coloured	10,0	1.100 x 1.100
010597020*	skin coloured	2,0	1.100 x 1.100
010597030*	skin coloured	3,0	1.100 x 1.100
010597040*	skin coloured	4,0	1.100 x 1.100
010597050*	skin coloured	5,0	1.100 x 1.100
010584020	red	2,0	1.100 x 1.100
010584030	red	3,0	1.100 x 1.100
010584040	red	4,0	1.100 x 1.100
010584050	red	5,0	1.100 x 1.100
010584070	red	7,0	1.100 x 1.100
010586020*	red	2,0	1.100 x 1.100
010586030*	red	3,0	1.100 x 1.100
010586040*	red	4,0	1.100 x 1.100
010586050*	red	5,0	1.100 x 1.100
010583050	yellow	5,0	1.100 x 1.100
010594010	black	1,0	1.100 x 1.100
010594020	black	2,0	1.100 x 1.100
010594030	black	3,0	1.100 x 1.100
010594040	black	4,0	1.100 x 1.100
010594050	black	5,0	1.100 x 1.100
010596020*	black	2,0	1.100 x 1.100
010596030*	black	3,0	1.100 x 1.100
010596040*	black	4,0	1.100 x 1.100
010591030	brown	3,0	1.100 x 1.100
010591040	brown	4,0	1.100 x 1.100
010593020*	brown	2,0	1.100 x 1.100
010593030*	brown	3,0	1.100 x 1.100
010593040*	brown	4,0	1.100 x 1.100
010593050*	brown	5,0	1.100 x 1.100

Further kind of Multiform with different Shore A and further colours available.

**Plastazote LD 45** (sheet goods)

Polyethylene cellular material, soft-elastic, closed-cell, light weight, very good for further processing, thermoforming at: approx. 160°C, hardness: approx. 15° Shore A, weight: approx. 0.045 g/cm<sup>3</sup>

**application:**

foot orthotics and footbeds, sport shoes



\* = perforated

Art. no.	Colour	Thickness approx. mm	Size approx. mm
010600020	white	2,0	1.000 x 1.000
010600030	white	3,0	1.000 x 1.000
010600040	white	4,0	1.000 x 1.000
010600050	white	5,0	1.000 x 1.000
010600060	white	6,0	1.000 x 1.000
010600100	white	10,0	1.000 x 1.000
010600200	white	20,0	1.000 x 1.000
010654020	black	2,0	1.000 x 1.000
010654030	black	3,0	1.000 x 1.000
010649030	grey	3,0	1.000 x 1.000
010649040	grey	4,0	1.000 x 1.000
010649050	grey	5,0	1.000 x 1.000
010649060	grey	6,0	1.000 x 1.000
010649100	grey	10,0	1.000 x 1.000
010605020	skin coloured	2,0	1.000 x 1.000
010605030	skin coloured	3,0	1.000 x 1.000
010605040	skin coloured	4,0	1.000 x 1.000
010605060	skin coloured	6,0	1.000 x 1.000
010630020*	skin coloured	2,0	1.000 x 1.000
010630030*	skin coloured	3,0	1.000 x 1.000

**HB-Foam** (sheet goods)

Polyethylene cellular material, very soft and light, very good for further processing

**application:**

foot orthotics and footbeds



Art. no.	Colour	Thickness approx. mm	Size approx. mm
010663203	grey	3,0	1.000 x 1.000
010663205	grey	5,0	1.000 x 1.000

**Savilux** (sheet goods)

EVA foam, extremely light weight, excellent thermoformable, easy handling, for making shoe insoles and bedding constructions, thermoforming at: approx. 90°–110°C, hardness: approx. 15° Shore A, weight: approx. 0.16 g/cm<sup>3</sup>

**application:**

foot orthotics and footbeds



Art. no.	Colour	Thickness approx. mm	Size approx. mm
010691030	skin coloured	3,0	1.000 x 560
010691040	skin coloured	4,0	1.000 x 560
010691060	skin coloured	6,0	1.000 x 560
010691080	skin coloured	8,0	1.000 x 560
010691100	skin coloured	10,0	1.000 x 560
010691120	skin coloured	12,0	1.000 x 560
010691140	skin coloured	14,0	1.000 x 560
010691200	skin coloured	20,0	1.000 x 560
010691240	skin coloured	24,0	1.000 x 560

Further colours and thicknesses available upon request.

# Schein padding and bedding materials

## Marilon (sheet goods)

Padding materials of EVA light cell, closed-cell, light weight, easily ground thermoforming at: approx. 100°–130°C, hardness: approx. 35°/40° Shore A, weight: approx. 0.16 g/qm<sup>3</sup>

### application:

foot orthotics and footbeds



Art. no.	Colour	Thickness approx. mm	Size approx. mm
010692040	green	4,0	600 x 940
010692060	green	6,0	600 x 940
010692080	green	8,0	600 x 940
010692100	green	10,0	600 x 940
010692120	green	12,0	600 x 940
010692140	green	14,0	600 x 940
010692200	green	20,0	600 x 940
010692240	green	24,0	600 x 940

## nora® Lunairmed (sheet goods)

EVA light cell, two-sided ginded, **specially suited for diabetics**, very soft, high elastically, good resilience, high walking comfort, closed-cell, washable, thermoforming at: approx. 120–170°C, hardness: approx. 18° Shore A, weight: approx. 0.08 g/cm<sup>3</sup>

### application:

footbeds, foot orthotics, forming components, padding plate for inner shoes, footbeds for rheumatics and diabetics



\* = perforated

Art. no.	Colour	Thickness approx. mm	Size approx. mm
011030503	skin coloured	3,0	1.050 x 825
011030504	skin coloured	4,0	1.050 x 825
011030505	skin coloured	5,0	1.050 x 825
011030506	skin coloured	6,0	1.050 x 825
011030508	skin coloured	8,0	1.050 x 825
011030510	skin coloured	10,0	1.050 x 825
011030512	skin coloured	12,0	1.050 x 825
011030524	skin coloured	24,0	1.050 x 825
011032903*	skin coloured	3,0	1.080 x 825
011032906*	skin coloured	6,0	1.080 x 825

## nora® Lunairflex (sheet goods)

EVA light cell, two-sided grinded, **specially suited for diabetics**, soft, high elastically, good resilience, closed-cell, washable, thermoforming at: approx. 120–170°C, hardness: approx. 22° Shore A, weight: approx. 0.12 g/cm<sup>3</sup>

### application:

footbeds, foot orthotics, padding plate for inner shoes etc.



\* = perforated

Art. no.	Colour	Thickness approx. mm	Size approx. mm
011035502	skin coloured	2,0	1.200 x 750
011035503	skin coloured	3,0	1.200 x 750
011035504	skin coloured	4,0	1.200 x 750
011035505	skin coloured	5,0	1.200 x 750
011035506	skin coloured	6,0	1.200 x 750
011035508	skin coloured	8,0	1.200 x 750
011035510	skin coloured	10,0	1.200 x 750
011035512	skin coloured	12,0	1.200 x 750
011036030*	skin coloured	3,0	1.200 x 750
011036060*	skin coloured	6,0	1.200 x 750

## nora® Lunalastic (sheet goods)

EVA light cell, two-sided grinded, excellent cushioning, high resilience, resistant to tearing, closed-cell, washable, thermoforming at: approx. 120–170°C, hardness: approx. 25° Shore A, weight: approx. 0.23 g/cm<sup>3</sup>

### application:

footbeds, foot orthotics, for soft-padded foot supports, padding plate for heel spur



\* = perforated

Art. no.	Colour	Thickness approx. mm	Size approx. mm
011050502	skin coloured	2,0	1.150 x 750
011050503	skin coloured	3,0	1.150 x 750
011050504	skin coloured	4,0	1.150 x 750
011050506	skin coloured	6,0	1.150 x 750
011050508	skin coloured	8,0	1.150 x 750

Further colours and thicknesses available upon request.

**nora® Lunasoft SLW** (sheet goods)

EVA light cell, ground both sides, flexible, elastic, washable and hygienic condition because of closed cellular structure,

thermoforming at: approx. 120°–170°C, hardness:

approx. 30° Shore A, weight: approx. 0.20 g/cm<sup>3</sup>

**application:**

footbed, foot orthotics and base construction

**nora® Lunasoft SL** (sheet goods)

EVA light cell, ground both sides, flexible, elastic, very high resilience, high walking comfort, washable and hygienic condition because of closed cellular structure,

thermoforming at: approx. 120°–170°C, hardness:

approx. 40° Shore A, weight: approx. 0.20 g/cm<sup>3</sup>

**application:**

footbed, foot orthotics and base construction



Art. no.	Colour	Thickness approx. mm	Size approx. mm
011045520	skin coloured	2,0	1.150 x 750
011045530	skin coloured	3,0	1.150 x 750
011045540	skin coloured	4,0	1.150 x 750
011045550	skin coloured	5,0	1.150 x 750
011045560	skin coloured	6,0	1.150 x 750
011045580	skin coloured	8,0	1.150 x 750
011045600	skin coloured	10,0	1.150 x 750
011045620	skin coloured	12,0	1.150 x 750
R10265303	white	3,0	900 x 650
R10265304	white	4,0	1.150 x 750
R10265306	white	6,0	1.150 x 750
R10265308	white	8,0	1.150 x 750
R10265310	white	10,0	1.150 x 750
R10265312	white	12,0	1.150 x 750
R10265314	white	14,0	890 x 640
R10265320	white	20,0	1.150 x 750
R10265408	med. brown	8,0	1.150 x 750
R10265410	med. brown	10,0	1.150 x 750
R10265508	black brown	8,0	890 x 640
R10265510	black brown	10,0	890 x 640
R10265512	black brown	12,0	890 x 640
R10265520	black brown	20,0	890 x 640
R10265602	black	2,0	1.150 x 750
R10265604	black	4,0	1.150 x 750
R10265606	black	6,0	1.150 x 750
R10265608	black	8,0	1.150 x 750
R10265610	black	10,0	1.150 x 750
R10265612	black	12,0	1.150 x 750
R10265620	black	20,0	1.150 x 750
R10265910	stone grey	10,0	1.150 x 750
R10266706	stone	6,0	1.150 x 750
R10266708	stone	8,0	1.150 x 750
R10266710	stone	10,0	1.150 x 750
R10266712	stone	12,0	1.150 x 750

Art. no.	Colour	Thickness approx. mm	Size approx. mm
011057020	skin coloured	2,0	1.140 x 740
011057030	skin coloured	3,0	1.140 x 740
011057040	skin coloured	4,0	1.140 x 740
011057050	skin coloured	5,0	1.140 x 740
011057060	skin coloured	6,0	1.140 x 740
011057080	skin coloured	8,0	1.140 x 740
011057100	skin coloured	10,0	1.140 x 740
011057120	skin coloured	12,0	1.140 x 740
R10269103	black brown	3,0	1.140 x 740
R10269106	black brown	6,0	1.140 x 740
R10269108	black brown	8,0	1.140 x 740
R10269110	black brown	10,0	1.140 x 740
R10269114	black brown	14,0	1.140 x 740
R10269120	black brown	20,0	1.140 x 740
R10269204	black	4,0	1.140 x 740
R10269206	black	6,0	1.140 x 740
R10269208	black	8,0	1.140 x 740
R10269210	black	10,0	1.140 x 740
R10269212	black	12,0	1.140 x 740
R10269214	black	14,0	1.140 x 740
R10269220	black	20,0	1.140 x 740
R10269224	black	24,0	1.140 x 740
R10269304	white	4,0	1.140 x 740
R10269306	white	6,0	1.140 x 740
R10269308	white	8,0	1.140 x 740
R10269310	white	10,0	1.140 x 740
R10269312	white	12,0	1.140 x 740
R10269314	white	14,0	1.140 x 740
R10269320	white	20,0	1.140 x 740
R10269406	med. brown	6,0	1.140 x 740
R10269408	med. brown	8,0	1.140 x 740
R10269410	med. brown	10,0	1.140 x 740
R10269420	med. brown	20,0	1.140 x 740
10270706	stone	6,0	1.140 x 740
10270708	stone	8,0	1.140 x 740
10270710	stone	10,0	1.140 x 740

Further colours and thicknesses available upon request.

# Schein padding and bedding materials

## ARU®batik 35 (sheet and roll goods)

Cover and padding EVA with batik effect, two sided-grinded, low specific gravity, increased tread damping with increased walking comfort, very attractive, thermoforming at: approx. 90°–140°C, hardness: approx. 35°–40° Shore A

### application:

footbeds, foot orthotics, specially for sport and children's shoes



\* = perforated

Art. no.	Colour	Thickness approx. mm	Size approx. mm
<b>sheet goods</b>			
R10337961	pink/white	1,0	1.100 x 950
R10337962	pink/white	2,0	1.100 x 950
R10337971	green/black	1,0	1.100 x 950
R10337972	green/black	2,0	1.100 x 950
R10337973	green/black	3,0	1.100 x 950
R10337991	black/red	1,0	1.100 x 950
R10337992	black/red	2,0	1.100 x 950
R10337993	black/red	3,0	1.100 x 950
R10338501	purple/green/black	1,0	1.150 x 875
R10338502	purple/green/black	2,0	1.150 x 875
R10339002	green	2,0	1.150 x 900
R10339003	green	3,0	1.150 x 900
R10339005	green	5,0	1.150 x 900
010639010	yellow/orange	1,0	1.100 x 950
010639020	yellow/orange	2,0	1.100 x 950
011639020*	yellow/orange	2,0	1.100 x 950
010639030	yellow/orange	3,0	1.100 x 950
010639040	yellow/orange	4,0	1.100 x 950
010639050	yellow/orange	5,0	1.100 x 950
010639100	yellow/orange	10,0	1.100 x 950
010641010	blue	1,0	1.030 x 1.035
010641020	blue	2,0	1.030 x 1.035
011641020*	blue	2,0	1.030 x 1.035
010641030	blue	3,0	1.030 x 1.035
010641040	blue	4,0	1.030 x 1.035
010641050	blue	5,0	1.030 x 1.035
010641060	blue	6,0	1.030 x 1.035
010641080	blue	8,0	1.030 x 1.035
010641100	blue	10,0	1.030 x 1.035
010643010	black	1,0	1.030 x 1.035
010643020	black	2,0	1.030 x 1.035
011643020*	black	2,0	1.030 x 1.035

Art. no.	Colour	Thickness approx. mm	Size approx. mm
<b>sheet goods</b>			
010643030	black	3,0	1.030 x 1.035
010643040	black	4,0	1.030 x 1.035
010643050	black	5,0	1.030 x 1.035
010643060	black	6,0	1.030 x 1.035
010643080	black	8,0	1.030 x 1.035
010643100	black	10,0	1.030 x 1.035
010638010	azure	1,0	1.030 x 1.035
010638020	azure	2,0	1.030 x 1.035
011638020*	azure	2,0	1.030 x 1.035
010638030	azure	3,0	1.030 x 1.035
010638040	azure	4,0	1.030 x 1.035
010638050	azure	5,0	1.030 x 1.035
010638060	azure	6,0	1.030 x 1.035
010638100	azure	10,0	1.030 x 1.035
010642010	beige	1,0	1.030 x 1.035
010642020	beige	2,0	1.030 x 1.035
011642020*	beige	2,0	1.030 x 1.035
010642030	beige	3,0	1.030 x 1.035
010642040	beige	4,0	1.030 x 1.035
010642050	beige	5,0	1.030 x 1.035
010642060	beige	6,0	1.030 x 1.035
010642100	beige	10,0	1.030 x 1.035

Art. no.	Colour	Thickness approx. mm	Size approx.
<b>roll goods</b>			
R10339801	orange	2,0	12 cm x 12 lfm.
R10339802	blue	2,0	12 cm x 12 lfm.
R10339804	black	2,0	12 cm x 12 lfm.
R10339806	beige	2,0	12 cm x 12 lfm.
R10339808	pink/white	2,0	12 cm x 12 lfm.
R10339810	green/black	2,0	12 cm x 12 lfm.
R10339814	black/red	2,0	12 cm x 12 lfm.

**Thermofix** (sheet goods)

EVA light cell, two sided-grinded, bodyfriendly, extremely thermoformable, light weight, stable in shape, easy so handle, thermoforming at: approx. 90°–100°C, hardness: approx. 45° Shore A

**application:**

footbeds, foot orthotics, bathing slippers, sportshoes

**Pryx** (sheet goods)

EVA light cell, two sided-grinded, bodyfriendly, extremely thermoformable, light weight, stable in shape, easy so handle, thermoforming at: approx. 90°–100°C, hardness: approx. 45° Shore A

**application:**

footbeds, foot orthotics, bathing slippers, sportshoes




Art. no.	Colour	Thickness approx. mm	Size approx. mm
010607302	white	2,0	1.560 x 900
010607303	white	3,0	1.560 x 900
010607304	white	4,0	1.560 x 900
010607306	white	6,0	1.560 x 900
010607308	white	8,0	1.560 x 900
010607310	white	10,0	1.560 x 900
010607312	white	12,0	1.560 x 900
010607320	white	20,0	1.560 x 900
010699303	brown	3,0	1.560 x 900
010699306	brown	6,0	1.560 x 900
010699308	brown	8,0	1.560 x 900
010699310	brown	10,0	1.560 x 900
010699314	brown	14,0	1.560 x 900
010699320	brown	20,0	1.560 x 900

Art. no.	Colour	Thickness approx. mm	Size approx. mm
011060030	white	3,0	0,5 m <sup>2</sup>
011060040	white	4,0	0,5 m <sup>2</sup>
011060060	white	6,0	0,5 m <sup>2</sup>
011060080	white	8,0	0,5 m <sup>2</sup>
011060100	white	10,0	0,5 m <sup>2</sup>
011060140	white	14,0	0,5 m <sup>2</sup>
011060200	white	20,0	0,5 m <sup>2</sup>
011062030	corky	3,0	0,5 m <sup>2</sup>
011062040	corky	4,0	0,5 m <sup>2</sup>
011062060	corky	6,0	0,5 m <sup>2</sup>
011062100	corky	10,0	0,5 m <sup>2</sup>
011062140	corky	14,0	0,5 m <sup>2</sup>
011062200	corky	20,0	0,5 m <sup>2</sup>



Further colours and thicknesses available upon request.

# Schein padding and bedding materials



## X2 soft medium (sqm goods)

PU foam, one-sided rough, delayed speed of returning to original position, excellent soft bedding, very low own weight, kind to the skin, easy to handle  
hardness: approx. 35° Shore A at 20°C

### application:

footbeds, foot orthotics, specially for diabetic treatment

\* = one side rough

\*\* = two-sided smooth

Art. no.	Colour	Thickness approx. mm	Size approx. mm
sqm goods			
011096030*	petrol	3,0	1.370 x 1.000
011096060**	petrol	6,0	1.370 x 1.000

## X2 soft (sqm and roll goods)

PU foam, one-sided rough, delayed speed of returning to original position, excellent soft bedding, very low own weight, kind to the skin, easy to handle

hardness: approx. 25° Shore A at 20°C

### application:

footbeds, foot orthotics, specially for diabetic treatment



\* = one side rough

\*\* = two-sided smooth

Art. no.	Colour	Thickness approx. mm	Size approx. mm
sqm goods			
011091030*	maroon	3,0	1.370 x 1.000
011091060**	maroon	6,0	1.370 x 1.000

Art. no.	Colour	Thickness approx. mm	Size approx. mm
roll goods			
R31109031*	maroon	3,0	20 m x 11 cm
R31109061**	maroon	6,0	10 m x 11 cm

## ARU®foam salmon (sqm and roll goods)

PU foam, soft, extra-light and airy, delayed speed of returning to original position, excellent soft bedding, kind to the skin, easy to handle, one-sided rough

### application:

footbeds, foot orthotics, specially for diabetic treatment



\* = one side rough

\*\* = two-sided smooth

Art. no.	Colour	Thickness approx. mm	Size approx. mm
sqm goods			
011121030*	salmon	3,0	1.370 x 1.000
011121045*	salmon	4,5	1.370 x 1.000
011121061*	salmon	6,0	1.370 x 1.000
011121060**	salmon	6,0	1.370 x 1.000

Art. no.	Colour	Thickness approx. mm	Size approx. mm
roll goods			
R31106031*	salmon	3,0	20 m x 11 cm
R31106061*	salmon	6,0	10 m x 11 cm



**P2 (sqm and roll goods)**

PU cellular material, highest ability to return to original position, air-permeable, can take permanent loads, kind to the skin, form stable, very good for further processing,

hardness: approx. 10°–12° Shore A,

**application:**

foot orthotics and footbeds, for many cushioning applications in orthopaedics and orthopaedic shoemaking



\* = two-sided rough

\*\* = perforated/one side rough

\*\*\* = one side rough

Art. no. sqm goods	Colour	Thickness approx. mm	Size approx. mm
011078015*	skin coloured	1,5	1.370 x 1.000
011072515**	skin coloured	1,5	1.370 x 1.000
011072015***	skin coloured	1,5	1.370 x 1.000
011078020*	skin coloured	2,0	1.370 x 1.000
011078030*	skin coloured	3,0	1.370 x 1.000
011072530**	skin coloured	3,0	1.370 x 1.000
011072030***	skin coloured	3,0	1.370 x 1.000
011078045*	skin coloured	4,5	1.370 x 1.000
011072045***	skin coloured	4,5	1.370 x 1.000
011078060*	skin coloured	6,0	1.370 x 1.000
011072060***	skin coloured	6,0	1.370 x 1.000
011072090***	skin coloured	9,0	1.370 x 1.000

Art. no. roll goods	Colour	Thickness approx. mm	Size approx. mm
R31107230***	skin coloured	3,0	20 m x 11 cm
R31107260***	skin coloured	6,0	10 m x 11 cm

**Scheinolon (sqm and roll goods)**

PU cellular material, very good to return to original position, excellent cushioning, form stable, very good for further processing,

hardness: approx. 10°–12° Shore A,

**application:**

foot orthotics and footbeds, for many cushioning applications in orthopaedics and orthopaedic shoemaking



\* = two-sided rough

\*\* = perforated/one side rough

\*\*\* = two-sided smooth

\*\*\*\* = one side rough

Art. no. sqm goods	Colour	Thickness approx. mm	Size approx. mm
011101030****	skin coloured	3,0	1.370 x 1.000
011103030*	skin coloured	3,0	1.370 x 1.000
011101530**	skin coloured	3,0	1.370 x 1.000
011101045****	skin coloured	4,5	1.370 x 1.000
011101060***	skin coloured	6,0	1.370 x 1.000
011101061****	skin coloured	6,0	1.370 x 1.000

Art. no. roll goods	Colour	Thickness approx. mm	Size approx. mm
R31110061****	skin coloured	6,0	10 m x 11 cm

**PX-Absorb (sheet goods)**

PU foam with delayed speed of returning to original position, soft, very durable, high shock absorbing

**application:**

footbeds, foot orthotics, specially for diabetic treatment



\* = one side rough

Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
011107030*	yellow	3,0	1.370 x 1.000
011107045*	yellow	4,5	1.370 x 1.000
011107060*	yellow	6,0	1.370 x 1.000

# Schein padding and bedding materials


## OnDry® by onSteam® embossed Poron/P2

(sqm and roll goods)

This material combination of Poron padding material and OnDry® by onSteam® micro fibre connects the excellent characteristics of these two materials. 100% breathable, regulating the temperature, odour resistant, easy to process, washable, best wearing comfort

### application:

footbeds, foot orthotics



Art. no.	Colour	Thickness approx. mm	Size approx. mm
R31105440	black	1,5	1.400 x 1.000
R31105443	black	3,0	1.400 x 1.000
R31105460	skin coloured	1,5	1.400 x 1.000
R31105463	skin coloured	3,0	1.400 x 1.000

Art. no. roll goods	Colour	Thickness approx. mm	Size approx. mm
R31105450	black	1,5	20 m x 11 cm
R31105453	black	3,0	20 m x 11 cm
R31105470	skin coloured	1,5	20 m x 11 cm
R31105473	skin coloured	3,0	20 m x 11 cm

## ARU® X-Static honeycomb with P2 (sqm goods)

Silver ions fabric laminated with P2, padding-, footbed and cover material, antibacterial, odour resistant, thermodynamical, antistatic, tear resistant, combines the features of P2 and X-Static, saves time and costs

### application:

footbeds, foot orthotics



Art. no. sqm goods	Colour	Thickness approx. mm	Size approx. mm
R10900288	black	1,5	1.400 x 1.000

## P2/PiGreco by schein

(sheet and roll goods)

Excellent material combination of P2 and the new microfiber "PiGreco by schein" to use for insoles and foot beddings. This padding, bedding and cover material combines the characteristics of P2 and PiGreco by schein and saves time and costs.



Art. no.	Colour	Thickness approx. mm	Size approx. mm
011083031	beige	3,0	1.370 x 1.100
019030020	black	2,0	1.370 x 1.100

Art. no. roll goods	Colour	Thickness approx. mm	Size approx. mm
R31105701	beige	1,0	20 m x 11 cm
R31105703	beige	2,0	20 m x 11 cm
R31105711	black	1,0	20 m x 11 cm
R31105713	black	2,0	20 m x 11 cm

## ARU® X-Static with P2 (sqm and roll goods)

Silver ions fabric laminated with P2, padding-, footbed and cover material, antibacterial, odour resistant, thermodynamical, antistatic, tear resistant, combines the features of P2 and X-Static, saves time and costs

### application:

footbeds, foot orthotics



Art. no. sqm goods	Colour	Thickness approx. mm	Size approx. mm
R10900290	black	1,5	1.500 x 1.000

Art. no. roll goods	Colour	Thickness approx. mm	Size approx. mm
R10900282	black	1,5	25 m x 11 cm

Further colours and thicknesses available upon request.

**ARU®foam salmon Dryspacer**

(sqm and roll goods)

PU foam with Dryspacer lamination, padding-, footbed and cover material, combines the features of ARU®foam salmon and Dryspacer, saves time and costs

**application:**

footbeds, foot orthotics,  
specially for diabetic treatment



Art. no.	Colour	Thickness approx. mm	Size approx. mm
sqm goods			
R31106135	black	3,0	1.370 x 1.000
R31106140	beige	3,0	1.370 x 1.000

Art. no.	Colour	Thickness approx. mm	Size approx. mm
roll goods			
R31106136	black	3,0	20 m x 11 cm

**ARU®foam salmon Spacer**

(sqm goods)

PU foam with Drylex lamination Spacer, padding-, footbed and cover material, combines the features of ARUfoam salmon and Drylex Spacer, saves time and costs

**application:**

footbeds, foot orthotics,  
specially for diabetic treatment



Art. no.	Colour	Thickness approx. mm	Size approx. mm
sqm goods			
R31106343	black	4,3	1.370 x 1.000

**P2/X2 Combi (sqm goods)**

PU foam combination, combines all of the features of P2 and X2, saves time and costs, hardness P2: approx. 10°–12° Shore A

**application:**

footbeds, foot orthotics,  
specially in the treatment of diabetes/rheumatism




Art. no.	Colour	Thickness approx. mm	Size approx. mm
sqm goods			
011076130	skin coloured 3,0 mm + maroon 3,0 mm	6,0	1.370 x 1.000

**P2/Poron+PE (Plastocom) (sheet goods)**

PU-Polyethylene combination, combines all features of P2/Poron and PE, saves time and costs hardness: approx. 14°–15° Shore A

**application:**

footbeds, foot orthotics



Art. no.	Colour	Thickness approx. mm	Size approx. mm
sheet goods			
010658250	blue PU 2,0 mm/ + skin coloured PE 3,0 mm	5,0	1000 x 1.000

# Schein/ARU combination plates

## ARUcombi DuoDirect®1 (sheet goods)

Padding plate, EVA light zell, directly vulcanized connection of 6 mm highly flexible, good resilience

**ARU®flexlastic**, approx. 25° Shore A,

and 10 mm very light, stabilizing


**ARU®soft AirLight**, approx. 50/55° Shore A,

thermoforming at: approx. 120°–130°C,

### application:

footbeds, foot orthotics, specially for diabetic treatment

6 mm **ARU®flexlastic**, ca. 25° Shore A  
10 mm **ARU®soft AirLight**, ca. 50/55° Shore A



**95 kg**

patients up to

Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340301	skin coloured/ white	16,0	950 x 570
10340372	skin coloured/ white	16,0	950 x 1.700

## ARUcombi DuoDirect®3 (sheet goods)

Padding plate, EVA light zell, directly vulcanized

connection of 6 mm highly flexible, good resilience

**ARU®flexlastic**, approx. 25° Shore A,

and 10 mm flexible, elastically


**ARU®flexis**, approx. 35° Shore A,

thermoforming at: approx. 120°–130°C,

### application:

footbeds, foot orthotics, specially for diabetic treatment

6 mm **ARU®flexlastic**, ca. 25° Shore A  
10 mm **ARU®flexis**, ca. 35° Shore A



**70 kg**

patients up to

Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340303	skin coloured/ grey	16,0	950 x 600
10340374	skin coloured/ grey	16,0	950 x 1.800

## ARUcombi DuoDirect®2 (sheet goods)

Padding plate, EVA light zell, directly vulcanized

connection of 6 mm soft padding

**ARU®liso**, approx. 20° Shore A,

and 10 mm flexible, elastically


**ARU®flexis**, approx. 35° Shore A,

thermoforming at: approx. 120°–130°C,

### application:

footbeds, foot orthotics, specially for diabetic treatment

6 mm **ARU®liso**, ca. 20° Shore A  
10 mm **ARU®flexis**, ca. 35° Shore A



**70 kg**

patients up to

Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340302	yellow/ grey	16,0	1.000 x 600
10340373	yellow/ grey	16,0	1.000 x 1.800

## ARUcombi DuoDirect®4 (sheet goods)

Padding plate, EVA light zell, directly vulcanized

connection of 6 mm highly flexible, good resilience

**ARU®flexlastic**, approx. 25° Shore A,

and 3 mm soft padding


**ARU®liso**, approx. 20° Shore A,

thermoforming at: approx. 120°–130°C,

### application:

footbeds, foot orthotics, specially for diabetic treatment

6 mm **ARU®flexlastic**, ca. 25° Shore A  
3 mm **ARU®liso**, ca. 20° Shore A



**50 kg**

patients up to

Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340304	skin coloured/ yellow	9,0	600 x 1.000
10340375	skin coloured/ yellow	9,0	1.000 x 1.800



All materials might be punched into individual shapes. Please ask our sales team for further information.

### ARUcombi DuoDirect®5 (sheet goods)

Padding plate, EVA light zell, directly vulcanized connection of soft padding 6 mm **ARU®liso**, approx. 25° Shore A, and 10 mm **ARU®liso**, approx. 30° Shore A, thermoforming at: approx. 120°–130°C,

**application:**

footbeds, foot orthotics, specially for diabetic treatment



patients up to



Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340305	light beige/ grey	16,0	600 x 1.000
10340376	light beige/ grey	16,0	1.000 x 1.800

### ARUcombi DuoDirect®7 (sheet goods)

Padding plate, EVA light zell, directly vulcanized connection of soft padding 7 mm **ARU®liso**, approx. 35° Shore A, and 10 mm **ARU®liso**, approx. 50/55° Shore A, thermoforming at: approx. 120°–130°C,

**application:**

footbeds, foot orthotics, specially for diabetic treatment



patients up to



Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340307	light green/ grey	17,0	950 x 750
10340378	light green/ grey	17,0	950 x 1.700

### ARUcombi DuoDirect®6 (sheet goods)

Padding plate, EVA light zell, directly vulcanized connection of soft padding 6 mm **ARU®liso**, approx. 25° Shore A, and 10 mm **ARU®liso**, approx. 50° Shore A, thermoforming at: approx. 120°–130°C,

**application:**

footbeds, foot orthotics, specially for diabetic treatment



patients up to



Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340306	beige/ dark brown	16,0	950 x 750
10340377	beige/ dark brown	16,0	950 x 1.700

### ARUcombi TricoDirect®2 (sheet goods)

Padding plate, EVA light zell, directly vulcanized connection of **EVA**, 3 mm, approx. 18° Shore A, skin, **ARU®foam Med„M“** blue 3 mm and **EVA**, 3 mm, approx. 25° Shore A, white. The result of this combination is an excellent thermoplastic deformable bedding material which is **ideally suited for diabetics**, thermoforming at: approx. 90°–120°C,

**application:**

footbeds, foot orthotics, specially for diabetic treatment



patients up to

Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340202	skin coloured/ blue/white	9,0	1.000 x 1.400

All materials might be punched into individual shapes. Please ask our sales team for further information.


# Schein/ARU combination plates

## ARUcombi TricoDirect®8 (sheet goods)

Padding plate, EVA light zell,  
directly vulcanized connection of soft padding  
5 mm **ARU®liso**, approx. 20° Shore A,  
5 mm **ARU®liso**, approx. 35° Shore A,  
10 mm **ARU®liso**, approx. 45° Shore A,  
thermoforming at: approx. 120°–130°C,

### application:

footbeds, foot orthotics, specially for diabetic treatment



patients up to


Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340308	light beige/red/ white	20,0	950 x 550
10340108	light beige/red/ white	20,0	950 x 825
10340379	light beige/red/ white	20,0	950 x 1.650

## ARUcombi TricoDirect®10 (sheet goods)

Padding plate, EVA light zell,  
directly vulcanized connection of soft padding  
5 mm **ARU®liso**, approx. 20° Shore A,  
5 mm **ARU®liso**, approx. 35° Shore A,  
10 mm **ARU®liso**, approx. 55/60° Shore A,  
thermoforming at: approx. 120°–130°C,

### application:

footbeds, foot orthotics, specially for diabetic treatment



patients up to


Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340310	beige/blue/ yellow	20,0	950 x 550
10340110	beige/blue/ yellow	20,0	950 x 825
10340381	beige/blue/ yellow	20,0	1.650 x 950

## ARUcombi TricoDirect®9 (sheet goods)

Padding plate, EVA light zell,  
directly vulcanized connection of soft padding  
6 mm **ARU®liso**, approx. 30° Shore A,  
6 mm **ARU®liso**, approx. 35° Shore A,  
8 mm **ARU®liso**, approx. 50° Shore A,  
thermoforming at: approx. 120°–130°C,

### application:

footbeds, foot orthotics, specially for diabetic treatment



patients up to


Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340309	beige/green/ dark brown	20,0	950 x 550
10340380	beige/green/ dark brown	20,0	1.650 x 950

## ARUcombi TricoDirect®11 (sheet goods)

Padding plate, EVA light zell,  
directly vulcanized connection of soft padding  
5 mm **ARU®liso**, approx. 20° Shore A,  
5 mm **ARU®liso**, approx. 30° Shore A,  
10 mm **ARU®liso**, approx. 40° Shore A,  
thermoforming at: approx. 120°–130°C,

### application:

footbeds, foot orthotics, specially for diabetic treatment



patients up to

Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340311	light beige/ grey/light green	20,0	950 x 550
10340382	light beige/ grey/light green	20,0	1.650 x 950


All materials might be punched into individual shapes. Please ask our sales team for further information.

**ARUcombi TricoDirect®12** (sheet goods)

Padding plate, EVA light zell,  
directly vulcanized connection of soft padding  
5 mm **ARU®Iiso**, approx. 20° Shore A,  
5 mm **ARU®Reboflex**, approx. 25° Shore A,  
12 mm **ARU®Iiso**, approx. 35° Shore A,  
thermoforming at: approx. 120°–130°C,

**application:**

footbeds, foot orthotics, specially for diabetic treatment



patients up to


Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340312	beige/ celeste/ dark blue	22,0	950 x 550
10340112	beige/ celeste/ dark blue	22,0	950 x 825
10340383	beige/ celeste/ dark blue	22,0	1.650 x 950

**ARUcombi TricoDirect®14** (sheet goods)

Padding plate, EVA light zell,  
directly vulcanized connection of soft padding  
4 mm **ARU®Iiso**, approx. 25° Shore A,  
4 mm **ARU®Iiso**, approx. 55° Shore A,  
12 mm **ARU®Iiso**, approx. 45° Shore A,  
thermoforming at: approx. 120°–130°C,

**application:**

footbeds, foot orthotics, specially for diabetic treatment



patients up to


Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340314	beige/hueso/ red	20,0	950 x 550
10340385	beige/hueso/ red	20,0	1.650 x 950

**ARUcombi TricoDirect®13** (sheet goods)

Padding plate, EVA light zell,  
directly vulcanized connection of soft padding  
4 mm **ARU®Iiso**, approx. 20° Shore A,  
4 mm **ARU®Iiso**, approx. 25° Shore A,  
12 mm **ARU®Iiso**, approx. 45° Shore A,  
thermoforming at: approx. 120°–130°C,

**application:**

footbeds, foot orthotics, specially for diabetic treatment



patients up to

Art. no. sheet goods	Colour	Thickness approx. mm	Size approx. mm
10340313	beige/light green/ brown	20,0	950 x 550
10340384	beige/light green/ brown	20,0	1.650 x 950

All materials might be punched into individual shapes. Please ask our sales team for further information.



# Schein thermoplastics and cork

## Erkoflex

rubber-like, transparent, soft, permanently elastic, can be ground and edged, thermoforming at: approx. 70°–80°C, manufacture of ankle- and peroneal paralysis caps, for interior shoe technique, test shoes, shaft samples



Art. no. sqm goods	Colour	Thickness approx. mm	Size approx. mm
010560015	transparent	1,5	1.000 x 1.000
010560020	transparent	2,0	1.000 x 1.000
010560030	transparent	3,0	1.000 x 1.000
010560040	transparent	4,0	1.000 x 1.000
010560050	transparent	5,0	1.000 x 1.000

## Erkolen

PE plastic, transparent, elastic, tear-resistant, can be vacuum-formed, thermo-forming at: approx. 120°C, reinforcing sole for elegant inlays, shell material for foaming of lasts



Art. no.	Colour	Thickness approx. mm	Size approx. mm
R57510106	transparent	0,6	1 qm

## ARU®-Superflex cork

pressed cork with PU bond, very light, flexible, easily grindable in both directions, slightly coarse texture, density: approx. 0.20–0.22 g/cm<sup>3</sup> for beddings and foot orthotics



Art. no.	Thickness approx. mm	Size approx. mm
R10500702	2,0	1.000 x 500
R10500703	3,0	1.000 x 500
R10500704	4,0	1.000 x 500
R10500705	5,0	1.000 x 500
R10500706	6,0	1.000 x 500
R10500708	8,0	1.000 x 500
R10500710	10,0	1.000 x 500
R10500720	20,0	1.000 x 500
R10500730	30,0	1.000 x 500

## ARU®-Multicork

EVA-corkmixture, thermoforming at: approx. 140–160°C, hardness: approx. 60° Shore A



Art. no.	Thickness approx. mm	Size approx. mm
010802702	2,0	950 x 925
010802704	4,0	950 x 925
010802706	6,0	950 x 925
010802708	8,0	950 x 925
010802710	10,0	950 x 925
010802712	12,0	950 x 925
010802715	15,0	950 x 925
010802720	20,0	950 x 925
010802725	25,0	950 x 925

## ARU®-Elastic cork

pressed cork with PU bond, very flexible, very easily grindable, slightly medium texture, for beddings and foot orthotics



Art. no.	Thickness approx. mm	Size approx. mm
R10500672	2,0	1.000 x 500
R10500673	3,0	1.000 x 500
R10500674	4,0	1.000 x 500
R10500675	5,0	1.000 x 500
R10500676	6,0	1.000 x 500
R10500678	8,0	1.000 x 500
R10500680	10,0	1.000 x 500

All materials might be punched into individual shapes. Please ask our sales team for further information.

**ARU®fant**

EVA-corkmixture, vacuum formable, elastic, light, thermoforming at: approx. 120°–140°C, hardness: approx. 50° Shore A, for beddings and foot orthotics, build-up material for contractions



Art. no.	Thickness approx. mm	Size approx. mm
010835020	2,0	940 x 535
010835030	3,0	940 x 535
010835040	4,0	940 x 535
010835060	6,0	940 x 535
010835080	8,0	940 x 535
010835100	10,0	940 x 535
010835150	15,0	940 x 535

**ARU®fant colour**

microcork, vacuum formable, elastic, light, thermoforming at: approx. 120°C, hardness: approx. 50°/55° Shore A, for beddings and foot orthotics



Art. no.	Colour	Thickness approx. mm	Size approx. mm
R01005602	stone	2,0	1.000 x 800
R01005603	stone	3,0	1.000 x 800
R01005604	stone	4,0	1.000 x 800
R01005605	stone	5,0	1.000 x 800
R01005606	stone	6,0	1.000 x 800
R01005607	stone	7,0	1.000 x 800
R01005608	stone	8,0	1.000 x 800
R01005610	stone	10,0	1.000 x 800
R01005612	stone	12,0	1.000 x 800
R01005615	stone	15,0	1.000 x 800
R01005620	stone	20,0	1.000 x 800
R01005630	stone	30,0	1.000 x 800

**ARU®fant soft**

microcork, vacuum formable, elastic, light, soft, thermoforming at: approx. 120°C, hardness: approx. 45° Shore A, for beddings and foot orthotics



Art. no.	Colour	Thickness approx. mm	Size approx. mm
R10500802	white	2,0	950 x 1.070
R10500803	white	3,0	950 x 1.070
R10500804	white	4,0	950 x 1.070
R10500806	white	6,0	950 x 1.070
R10500808	white	8,0	950 x 1.070
R10500810	white	10,0	950 x 1.070
R10500812	white	12,0	950 x 1.070
R10500815	white	15,0	950 x 1.070

**Flexo cork**

flexible, pressed cork with conventional bonding, light, particularly elastic, easy to cut and grind, density: approx. 0.30–0.40 g/cm<sup>3</sup> for beddings, foot orthotics and last construction



Art. no.	Thickness approx. mm	Size approx. mm
010806020	2,0	1.000 x 500
010806030	3,0	1.000 x 500
010806040	4,0	1.000 x 500
010806050	5,0	1.000 x 500
010806060	6,0	1.000 x 500
010806080	8,0	1.000 x 500
010806100	10,0	1.000 x 500

**Superlight cork**

pressed cork with high-elasticity synthetic resin bonding, light, durable, bends in every direction, does not dry out even after extended storage density: approx. 0.25–0.35 g/cm<sup>3</sup> for beddings, foot orthotics, built-up wedges in shoe preparation



Art. no.	Thickness approx. mm	Size approx. mm
010860020	2,0	1.000 x 500
010860030	3,0	1.000 x 500
010860040	4,0	1.000 x 500
010860050	5,0	1.000 x 500
010860060	6,0	1.000 x 500
010860080	8,0	1.000 x 500
010860100	10,0	1.000 x 500
010860120	12,0	1.000 x 500
010860150	15,0	1.000 x 500

**Colasto cork**

pressed cork with high-elasticity synthetic resin bonding, light, durable, bends in every direction, does not dryout even after extended storage density: approx. 0.30–0.40 g/cm<sup>3</sup> for beddings, foot orthotics and last construction



Art. no.	Thickness approx. mm	Size approx. mm
010804010	1,0	1.000 x 500
010804020	2,0	1.000 x 500
010804030	3,0	1.000 x 500
010804040	4,0	1.000 x 500
010804050	5,0	1.000 x 500
010804060	6,0	1.000 x 500
010804080	8,0	1.000 x 500
010804100	10,0	1.000 x 500

Schein and its daughter company Ruckgaber offer the whole range of materials/tools/machinery/diagnostic devices for shoemaker's demand and for manufacturing prosthetics and orthotics.

For further information please follow: <https://www.schein.de/en/home/home.html>






COMMUNICATE WITH US  
VIA WHATSAPP:



**+49 2191 910-200**

[www.schein.de/whatsapp/privacy](http://www.schein.de/whatsapp/privacy)

04.0923.01.04-21



Schein Orthopädie Service KG  
Hildegardstr. 5  
42897 Remscheid, Germany  
Tel. +49 2191 910-200  
Fax +49 2191 910-100  
[service@schein.de](mailto:service@schein.de)  
[www.mygeneration-schein.de](http://www.mygeneration-schein.de)



**WeXhelp**

099993507

We reserve the right to make minor changes to colours and models.